

# AID 500 Series Drying case

For Rechargeable Hearing Aids

1 AID 500

2 AID 500-UV




## Feature

- This drying case is especial designed with lower temperature heating to meet the demand to dry rechargeable hearing instruments, provides an effective and affordable way to remove moisture.
- Heating temperature  $\leq 35^{\circ}\text{C}$
- Optional timer- 3/6 hours dry/ Touch key to switch
- Autostop- 30min UV
- Low voltage input - 5V
- Strong compatibility of power supply-2.0 USB cable
- Convenient to carry

## Six Reasons To Choose



## Suitable For Hearing Devices



# NEW UPGRADED DYING CASE

Touch Design    Drying&UV Sterilizing

