

- Used in combination with plug-in earphones.
- Comfortable to wear and easy to cause swelling.
- PU material, soft and delicate, perfect fit your ear canal.



- The foam expands after being placed in the ear canal to ensure a perfect fit.
- Fix the plug-in earphone during the whole measurement process to help get accurate hearing test results.



Delicate and soft material

Moderate rebound time, comfortable feeling to wear.

PU Material

Non-toxic and harmless environmental protection material.

Certification Assurance

REACH, RoHS, Test Report.



Safe material



Comfortable feeling to wear



Authority certification



Material: PU foam



Size: **Large**
Product size: 26mm(H)
Foam size: 13(H)*18(φ)mm



Size: **Middle**
Product size: 26mm(H)
Foam size: 13(H)*13(φ)mm



Size: **Small**
Product size: 26mm(H)
Foam size: 13(H)*10(φ)mm

INSERT



FOAM

TIP