

Catalogue-8



SOUNDLINK
Soundlink Co., Ltd.



2F-1 # Building North, 89 Songshan Road,
Shishan Industrial Complex Park, New District, Suzhou, China 215129

+86 512-69372012 info@soundlink.com.cn
+86 512-69008100-810 www.soundlink.cn

SOUNDLINK
SOUNDLINK AND MORE

ABOUT SOUNDLINK

Soundlink, established in 2004, is a global supplier for hearing aid accessories and solutions. It is located in Suzhou city where the cradle of hearing aid is in China.

With 19 years of rapid development, Soundlink has become a strategic partner of many domestic and international hearing industry providers, serviced to nearly 10,000 professional hearing organizations which from over 100 countries and areas.

We are focused on the field of hearing aid related solutions, noise protection and free hearing. With the original technology research & development and global resource integration, we are trying our best to meet the diverse needs of human hearing. Soundlink always strives to become a continuous leader in the development of the hearing industry!

Our corporate culture is "sincerity, seriousness, honesty, responsibility and gratitude". With rapid response, active service, reliable guarantee, Soundlink hopes to create a bright future of hearing business with you!

PRODUCT LINE

Hygiene & Care	01-21
Ear Impression	22-34
Marketing & Display	35-43
Earmold lab material & Equipment	44-62
Drilling & Grinding	63-69
Hearing aid accessories	70-82
Customized Earmolds	83-93
3D Printing Solutions	94-105



HYGIENE & CARE

Complete drying & cleaning product supplying for hearing instruments daily maintain.

Drying case

With patented hearing technology, daily removing moisture inside and on the surface of earmolds and hearing aids

Dimension: 80x80x80 (H)mm
Capacity: 70x70x19 mm
Package: 60pcs/carton
Temperature: 40°C ±5°C
Power: ≤ 8W
Material: ABS

CE Certificate



Aid 200

0201017-001 (Flat Plug)
0201017-002 (Round Plug)

Model: Aid200, no timer
Voltage: 220V



Aid 200T

0201018-001 (Flat Plug)
0201018-002 (Round Plug)

Model: Aid200T, 3h/6h timer
Voltage: 220V




AID 201

NEW



AID 201



AID 201-UV

Drying Case 201 Series

0201050-003 AID 201

0201050-004 AID 201-UV

Specially desinged for hearing aids drying.

Feature

- 1.Constant heating - 50 C ±5
- 2.Autostop - 8 hours dry, 30min UV
- 3.Low voltage input - 5V
- 4.Strong compatibility of power supply-2.0 USB cable
- 5.Easy to use
- 6.Convenient to carry
- 7.Size: 76.5 x 67 x 44mm
- 8.2 options -UV or without UV

FOUR REASONS TO CHOOSE



AID 202

NEW



Battery position Multipurpose brush

Drying Case 202 Series

0201040-002 (EU Plug)

0201040-003 (US Plug)

0201040-004 (UK Plug)

Multifunctional dryer for hearing aids.

Optional timer control: 3/6 hours


Temperature control: 40 °C /50 °C

(Users can switch freely according to actual needs)


Two positions

Battery position

Multipurpose brush



EIGHT REASONS TO CHOOSE



SUITABLE FOR HEARING DEVICES



AID 300A

NEW

Sanitizer Dryer UV LED AID 300A

- 0201039-001 (US Plug)
- 0201039-002 (EU Plug)
- 0201039-005 (UK Plug)



Introduction

- Upgraded from UV lamp to UV LED, with the latest design, and multiple functions, Aid 300A dryer is highly welcomed for hearing aids daily drying and disinfection.
- Battery test function: easy to read battery power.
- UV LED: longer life, 275nm UV C light for sterilization and disinfection.
- USB connector, reliable adaptor.

Technical Specifications

1.Certification:CE/SGS/UV Safety Report

Material:ABS
Outer Size:166x105x48mm
Inner Size:87x95mm
Net Weight:460g

2.Power

Voltage:DC 5.0V
Current:DC 2.0A
Power:≤5W

3.Automatic temperature control:50(±5 C)

4.Timer control


Optional drying time: 3h/6h /9h
Sterilize time:30 mins

5.Battery indicator:three levels

Green light-battery power ≥ 1.3V
Yellow light-battery power ≥ 1.15V but < 1.3V
Red light -battery power < 1.15V

6.Safety protection:

Over-temperature & over-heating protection
Protection against UV leaking



ZB03-B Vacuum Pump

- 0204004-001 (Flat Plug)
- 0204004-002 (Round Plug)
- 0204004-005 (Flat Plug/AC110V)

Introduction

- The device adopts the vacuum with dehumidifying drying technology which can dry thoroughly, and extend the service life of hearing aids.
- Air exhaust makes the negative pressure to absorb out the cerumen and dirt from cavity of the receiver and microphone. Airflow of the blowing air can help us clean the dust and foreign object of hearing aid surface.

Technical Specifications

1. Environment:
Ambient temperature: 5℃ - 40℃
Air humidity: ≤ 80%RH (when temperature is 25℃)
Atmospheric pressure: 86kpa~106kpa
No flammable, corrosive gas or conductive dust exist in the air
2. Power:
Voltage: AC 220V±10% / AC110±10%
Mains frequency: 50Hz±1Hz
Input: ≤ 100VA
Fuse: RFD5x20 1.5A/250V
3. Limited negative pressure: ≥-0.75Mpa

4. Vacuum drying tank dimension:Φ100x90(mm)
5. Timing accuracy:
Electrical:
The tolerance of test accuracy in sec. :±1s
The tolerance of test accuracy in min. :±3s
6. Noise: ≤ 65dB(A)
7. External dimension: 280x186x196 (H)mm
8. Weight: 5kg





ZB03-B II Vacuum Pump

0204005-001 (Flat Plug)

0204005-004 (Round Plug)




Feature

- Exquisite design, with user-friendly interface
- Convenient button operation, one-button smart starting
- Dehumidification maintenance, you only need this vacuum pump

Specification

- Time alarm: 15mins
- Working noise: ≤65dB
- Out size: 345 x 200 x 255(mm)
- Drying tank size: φ100x90(mm)
- Work power: 50Hz±1Hz
- Overload protection: Rfφ5x20 1.5A/250V

EXQUISITE DETAILS



Drying care instrument

NEW



Vacuum Desiccator EW01-M

0201014-002 (Round Plug)

0201014-003 (Flat Plug)



Feature

- LED sterilization lamp
Three LED lamp to sterilize all-around.
- Manual pressure relief valve
Double insurance
- Working surface
One-click to work
- Power switch
One button to turn off the power
- High transmittance ultraviolet protection cover
Protect the user's eyes and skin from ultraviolet damage observe in peace
- Certificate
CE
Safety Report of LED Sterilization Lamp
Effect Report of Sterilization and Disinfection

Date Sheet

Voltage:	DC 12V/5A	Relative air humidity:	≤80% RH
Working power:	≤40W	Extreme negative pressure:	≤-75kPa
Timing:	15 min	Temperature:	25 C
Working noise:	≤55dB	Atmospheric pressure:	86kPa-106kPa





0201006-006

Drying pallet - I

Used with drying jar for drying hearing aids and earmolds
The color will turn from orange to transparent when it is saturated with moisture.

Net weight: 4.5g
Dimension: Φ35x10 mm
Indicator: Silicone dioxide
Package: aluminum foil/pc ,20pcs/pack

CE certificate, RoHS certificate



0201022-002

Drying pallet - II

Used with drying jar for drying hearing aids and earmolds
The color will turn from orange to transparent when it is saturated with moisture.

Net weight: 30g
Dimension: Φ63x20mm
Indicator: Silicone dioxide

CE certificate, RoHS certificate



0201022-009

Drying pallet - II with box

Used with drying jar for drying hearing aids and earmolds
The color will turn from orange to transparent when it is saturated with moisture.

Package: 5pcs/box
Dimension: 105x68x68mm
Indicator: Silicone dioxide

CE certificate, RoHS certificate



0201035-001

Drying pallet - III

Used with drying jar for drying hearing aids and earmolds
The color will turn from orange to transparent when it is saturated with moisture.

Net weight:4.5g
Dimension: Φ31x12 mm
Indicator: Silicone dioxide
Package:2pcs/card, 50cards/carton

CE certificate, RoHS certificate



0201035-002

Drying pallet - III with box

Used with drying jar for drying hearing aids and earmolds
The color will turn from orange to transparent when it is saturated with moisture.

Package: 2cards/box
Indicator: Silicone dioxide

CE certificate, RoHS certificate



0201051-001

Reusable drying pallet

Used with drying jar for drying hearing aids and earmolds.
The color will turn from orange to green when it is saturated with moisture.
Can be used with microwave.

Dimension: φ62.5x20mm

CE certificate, RoHS certificate

Feature



Orange particles

- Irregular blocky structure
- Faster moisture absorption & Environmental protection



Discoloration

- Replace the drying pallet when the color changes to transparent.

Note




- Tear open the sealing aluminum foil when using
- Replace the drying pallet when the color changes to transparent.



- Do not tear off the hygroscopic paper.
- Do not eat .

Used with different type of drying jar

- Sequence of product placement
Drying pallet — Sponge — Hearing aid
- Close the lid tightly



Drying jar - I

Used with drying pallet(0201035-001) for drying hearing aids and earmolds
Plus sieve with slim vent

Material: PP
Dimension: Φ60x85 mm
Package: 50pcs/carton



Drying jar - II

Used with drying pallet for drying hearing aids and earmolds.
Especially designed for our drying pallet (0201022-002) to offer dehumidifier service for hearing instruments

Material: PP
Dimension: Φ72x68 mm
Package:48pcs/carton

Certificate: Biocompatibility report



Drying jar - III

Used with drying pallet(0201022-002) for drying hearing aids and earmolds.

Material: PP
Dimension: Φ68x54 mm
Package:48pcs/carton



Cleaning jar

Used with cleaning tablets for cleaning earmolds
Plus sieve with slim vent

Material: PP
Dimension: Φ60x85 mm
Package: 50pcs/carton



Cleaning tablets

Used with cleaning jar (0201011) for cleaning earmolds and killing bacteria and germs

Dimension: 138x92x5 (H)mm
Package: 60cards/carton



Cleaning tablets with box

Used with cleaning jar (0201011) for cleaning earmolds and killing bacteria and germs

Package: 6 cards / box



Cleaning wipes

Wipe and cleaning the surface of hearing aid & earmold.

Package: 100pcs/box
Composition: 70% alcohol



0205017-007

Cleaning wipes

Wipe and cleaning the surface of hearing aid & earmold.

Brand: Soundlink
Package: 30pcs/box
Composition: Aqua, Chlorhexidine diacetate, Cetrimonium bromide



0205018

Cleaning spray

Cleaning solution for earmold and shell of hearing aids

Brand: Nice-fit
Capacity: 30ml
Expiry date: 36 months
Package: 50bottles/carton
Contents: Cationic surfactants



0205001

Cleaning cloth

Soft, efficient to wipe earmolds and hearing aids
Color of cloth will be shipped randomly

Dimension: 140x145mm
Material: 80% terylene+20% chinlon



0205030

Cleaning cloth pocket

Soft and efficient to wipe and store earmolds and hearing aids,
to avoid friction and collision

Dimension: 105x90 mm
Material: 80%terylene + 20% chinlon




0205002

Cerumen stick

Ideal to remove moisture & earwax in hearing aids

Dimension: Φ35x195 mm

Do not use it to human ear



0202001

Brush with magnet

Used to clean hearing aids and take battery out with the magnet,
especially for the elderly

Package: 200pcs/pack
Handle Material: ABS
Brush Material: Nylon

Certificate: RoHS



0202002

Brush with screwdriver

Used to clean hearing aids and adjust trimmers

Package: 200pcs/pack
Handle Material: ABS
Brush Material: Nylon

Certificate: RoHS



0202003

Brush with loop&magnet- I

Used to clean hearing aids, take out battery and remove earwax in
sound outlet

Length: 50mm
5bunches
Package: 200pcs/pack
Handle Material: ABS
Brush Material: Nylon

Certificate: RoHS



0202008

Brush with loop&magnet- II

Used to clean hearing aids, take out battery and remove earwax in
sound outlet

Length: 61.5mm
7bunches
Package: 200pcs/pack
Handle Material: ABS
Brush Material: Nylon

Certificate: RoHS



0202009

Cleaning brush

Used to clean hearing aids and earmolds

Package: 200pcs/pack

0202011



Brush with wax loop

Used to clean wax or dirty on hearing aids and earmolds

5bunches
Package: 200pcs/pack
Handle Material: ABS
Brush Material: Nylon

Certificate: RoHS

0202005



Wax loop

Used to remove earwax or dirt in sound outlet

Package: 200pcs/pack
Handle Material: ABS

Certificate: RoHS

0202004



Vent cleaner

Used to clean vent of ITE,ITC,CIC

Package:500pcs/pack

Certificate: RoHS

0202014-001



Vent cleaner for open fitting

Used to clean vent of ITE,ITC,CIC and sound tube of open fitting hearing aid

Length: 12cm

0202014-002



Brush with magnet & thin spiral

Used to clean hearing aids and earmolds
With magnet on one side to take battery out and thin spiral on the other side

5bunches

0202014-003



Brush with magnet & thick spiral

Used to clean hearing aid and earmold.
With magnet on one side to take battery out and thick spiral on the other side

5bunches

0202014-004



Brush with loop, magnet & vent cleaner

Used to clean hearing aids and earmolds
With magnet on the top to take battery out vent cleaner to clean sound tube wax loop to remove earwax or dirt in sound outlet

9bunches

0202014-005



Screwdriver

Used to adjust trimmers

0202014-006



Remover for battery door

Used to remove the hearing aid battery door

0202014-007



Removable multipurpose brush

Removable cleaning brush with magnet, screwdriver and spiral.
Used to cleaning hearing aid, adjust trimmers and remove battery door

4bunches

Brush with screwdriver-II



0202014-008

Used to clean hearing aid and adjust trimmers

Brush with magnet & screwdriver




0202014-009

Used to clean hearing aids and adjust trimmers

5bunches

Air blower




0203001
0203001-001

Ideal to remove earwax and moisture in earmolds and tubing
Plastic opening tip

Material: PVC
Dimension: 40x100 mm
Package: 100pcs/carton

CE, RoHS, REACH

Air blower



0203003
0203003-001
0203003-002

Ideal to remove earwax and moisture in earmolds and tubing
Metal opening tip

Material: PVC
Dimension: 40x100 mm
Package: 100pcs/carton

CE, RoHS, REACH

0204002-002
0204002-005



Ultrasonic cleaner- I

Ideal to clean earmold, especially for hearing aid centers service

Power: 35W
Capacity: 600ml
Voltage: 0204002-002 Round plug(220V)
0204002-005 Flat plug(110V)
Frequency: 40,000 energy waves per second

0204003-002
0204003-016



Ultrasonic cleaner- II

Digital type, ideal partner to clean earmolds
Special and necessary for hearing aid centers service

Power: 50W
Capacity: 600ml
Frequency: 40KHz
Time control: 180s/360s
Material: stainless steel
Power supply: AC220-240V 50~60Hz
Type: 0204003-002 Flat plug
0204003-016 Round plug

0204003-003
0204003-017



Ultrasonic cleaner- III

Ultrasonic and heating design, ideal partner to clean earmolds
Special and necessary for hearing aid centers service

Capacity: 2L
Timer: 1-30mins
Frequency: 40KHz
Heating Power: 60W
Ultrasonic Power: 70W
Material: stainless steel
Power supply: AC220-240V 50~60Hz
Type: 0204003-003 Flat plug
0204003-017 Round plug

0204002-006



Ultrasonic cleaner-IV

With higher ultrasonic power conversion efficiency and faster cleaning speed, this ultrasonic cleaner makes your hearing aid center more meticulous, more intimate, and more professional

Brand: Soundlink
Capacity: 9L
Frequency: 40KHz
Heating Power:300W



Hearing Aid Clip for Children

Product Characteristic

- This product has up-to-date design and delicate workmanship with high quality.
- Twelve vivid and cute animals pattern on the plate of clip is attractive to children.
- The elastic rope and strong clip with humanized design protect BTE against falling off.
- Silicone adaptors are equipped with Biocompatibility Report.
- Good partner for Kids BTE.

Product Details



BTE clip (Binaural)

Prevent BTE from falling off, especially for children and seniors

Components: 1 x clamp, 1 x rope, 6 x adaptor

Length: 350mm

Adaptor dia: 7mm, 9mm, 11mm

Colour: 0205008-001	Grey	0205008-002	Brown
0205008-003	Dark blue	0205008-004	Light blue
0205008-005	Green	0205008-006	Black



BTE clip (Monaural)

Prevent BTE from falling off, especially for children and seniors

Components: 1 x clamp, 1 x rope, 3 x adaptor

Length: 350mm

Adaptor dia: 7mm, 9mm, 11mm

Colour: 0205006-001	Grey	0205006-002	Brown
---------------------	------	-------------	-------

BTE clip adaptor

Accessories of BTE clip

Material: Silicone

Dimension: F0205028-001	S	7mm
F0205009	M	9mm
F0205010	L	11mm



BTE jacket

Protect BTE against sweat and falling off

Material: Silicone


Model: 13#	0205019
675#	0205020

BTE jacket adaptor

Accessories of BTE jacket

Material: Silicone

Model: 13#	F0205013
675#	F0205014





0205009-012

Customized hearing aids BTE clip

Prevent customized hearing aids from falling off;
Elastic string to avoid pulling down and hurting

Color: Beige
Length of string: 360mm

Crocodile clip



0205029-001

Special accessories for BTE lanyard, suitable for clipping on thin collar for fixing the lanyard

Color: 0205029-001 Transparent

Hearing aid protector



0205024

Protect hearing aids against sweat and falling off
Random color for shipping

Material: Cotton



Drying kit



0205021-010

You can make up your own kit by adding any small items you prefer.

Components:
1x Drying jar
1x Air blower
1x Cleaning cloth
4x Drying pallet
1x Digital battery tester
1x Brush with loop&magnet- I

Cleaning kit



0205022-004

You can make up your own kit by adding any small items you prefer.

Components:
1x Cleaning jar
1x Air blower
1x Cleaning cloth
6x Cleaning tablets
1x Digital battery tester
1x Brush with loop&magnet-I



0104015

Clear kit bag- I

Used for carrying care-kit

Material: PVC
Dimension: 170x105x50 (H)mm



0104039-001

Clear kit bag-III

Used for carrying care-kit

Material: PVC
Color: Grey
Dimension: 155x65x135 (H)mm



Kit bag- I

Used for carrying care-kit

Material: Nylon oxford
Dimension: 165x65x135 (H)mm
Color: Grey (0104026-001)
Orange (0104026-002)



Kit bag- II

Used for carrying care-kit

Material: Nylon oxford
Dimension: 250x110x170 (H)mm
Color: Blue (0104025-001)
Green (0104025-002)

EAR IMPRESSION

Complete and necessary material
&
tools for impression taking and ear
canal checking



0306002-009

Impression-taking kit

Used to take ear impression
Case size: 330×280×80 (H)mm

Contents:

- 1.Otoscope- II
- 2.Impression syring- I
- 3.Impression foam: 2 packs (M size)
- 4.Tweezers: Bayonet shape
- 5.Ear light (LED) with probe tip and ear picker
- 6.Impression material: 400g×2 two jars with dosing spoons



0306003

Case of impression taking kit

Storage for impression taking tools

Material: Plastic

Dimension: 330×280×80 (H)mm



0401004

ABR impression material 700g×2

Mixing rate 1:1

With a very soft,slightly pressure-forming consistency,
high stability,impression exactitude,tensile strength.

Work well for BTE,ITE,ITC and CIC impressions.

Brand: ABR

Type: 700g×2

Color: white and green

Hardness: 37 shore±2



0401005

ABR impression material 400g×2

Mixing rate 1:1

With a very soft,slightly pressure-forming consistency
high stability,impression exactitude,tensile strength.

Work well for BTE,ITE,ITC and CIC impressions.

Brand: ABR

Type: 400g×2

Color: white and green

Hardness: 37 shore±2



0401017-001

ABR impression material 800g×2

Mixing rate 1:1
With a very soft,slightly pressure-forming consistency,
high stability,impression exactitude,tensile strength.
Work well for BTE,ITE,ITC and CIC impressions.

Brand: ABR Type: 800g×2
Color: white and green Scent: lemon
Hardness: 36 shore±2



0301002-002
0301002-003

Impression syringe - I

Reusable to inject impression material into ears

Material: PC Dimension: 3.5mm
Color:Blue-0301002-002 White-0301002-003



0401019

Green eco impression material 500g×2

High elastic recovery due to an innovative polymer,exceptional
tear strength,lasting dimensional stability, impression exactitude
Work well for BTE,ITE,ITC and CIC impressions

Brand: Detax
Type: 500g×2
Scent: lemon
Color: green and yellow
Hardness: 40 shore±2



0301013

Electric injector

Chargeable electrically and effortlessly operated injector
for dispensing the impression material from double artridges.
The acceleration speed and the dispensed amount of material stay
consistent.
It can also be adjusted to suit situational needs.

Color: white
Votage: 110V-240V
Mixing rate: 1:1



0401025


Dreve ear impression material 800g×2

High elastic recovery due to an innovative polymer,exceptional
tear strength,lasting dimensional stability, impression exactitude
Work well for BTE,ITE,ITC and CIC impressions

Brand: Dreve
Type: 800g×2
Color: green and white
Hardness: 35 shore±2



0301004-005
0301004-006



Manual inject pistol

For inserting standard double cartridges for impression taking.
By using the inject pistol,the impression material can be delivered to the
ear while applying the steady pressure.

Mixing rate: 1:1



0301022-004
0301023-002

Impression syringe - II

Reusable to inject impression material into ears

Material: PC
Dimension: 2.5 mm (0301022-004)
3.5 mm (0301023-002)



0301005

Mixing tip

For standard double cartridges for automatic mixing
of the impression material

Material: PP
Dimension: Φ7×78 mm
Package: 100pcs/bag



0301005-003

Mixing tip(MB 1:1)

For standard double cartridges for automatic mixing of the impression material

Material: PP
Dimension: Ø5.4 mm
Package: 100pcs/bag



0301006

Canula tip - I

Mixing canula to facilitate the impression taking of very narrow auditory canals. Suitable for the standard mixing canula

Material: PP
Package: 100pcs/bag



0301020

Canula tip - II

Mixing canula to facilitate the impression taking of very narrow auditory canals. Suitable for the standard mixing canula

Material: PP
Package: 100pcs/bag



0302033
0302034

Disposable specula

Works with otoscope(0302004-002)

Type: 2.5mm-0302033
4.0mm-0302034



0302003-007

Spare bulb for otoscope 0302004-002

Bulb type: LED



0302008-006

Ear light

With LED for easy viewing of the ear canal

Power supply: AAAx2
Accessory: probe tip, ear picker

Battery not included



F0302003
F0302002

Ear Light parts

Accessory for ear light(0302008-006)

Type: Ear Light Probe tip Length: 51mm-F0302003
Ear Light Picker Length: 53mm-F0302002



0302028

Ear light Riester

Color: Black 0302009-001
White 0302009-002



0302004-002

Otoscope

With optical fiber light for an excellent view of ear canal

Speculum Tip Size : Dia.2.5mm; 3mm; 4mm; 5mm
Lens Diameter : 18mm
Illumination : Direct Illumination
Luminance : 24 Lumen
LED Bulb : 108G 2.5V 750mAh
Color Temperature : 3050K +/-5%
Working Situation : Temperature : 10-35 °C; Humidity : 30-75%;
Air Pressure : 700-1060hPa

Battery not included

e-scope® Otoscopes

Product introduction

- Riester e-scope® otoscopes are available with halogen or vacuum direct illumination, F.O. xenon or F.O. LED
- Sealing precision optics pivotable to both sides with triple magnification for optimum view
- Sturdy ear speculum receptacle in hardwearing hygienic metal for e-scope® /ri-scope® L1/L2 speculum
- Compatible with ear speculae from other renowned manufacturers
- Suitable for performance of pneumatic tests (supplied without connector and insufflator)
- Available in a choice of black or white finish
- Economically priced fibre optic version with 2.5v Xenon (3200 K)*
- Fibre optics for optimum focusing and guiding of the light

Brand: Riester
 Voltage: LED: 3.7V
 Xenon: 2.5V
 Halogen: 2.5V
 Power supply : AAx2 batteries

Kelvin data: LED: 5500 K
 Xenon: 3200 K
 Halogen: 3000 K



e-scope® Otoscopes



□ White



■ Black

Halogen Bulb 2.5V Otoscope

Color : White 0302002-002
 Black 0302002-003



□ White



■ Black

Fibre Optic 2.5V Xenon Otoscope

Color : White 0302002-004
 Black 0302002-005



□ White



■ Black

LED Fibre Optic 3.7V Otoscope

Color : White 0302002-006
 Black 0302002-007

Light bulb



Halogen bulb
0302003-003



Xenon bulb
0302003-004






LED bulb
0302003-005



0302039-001

Miis Horus OTOSCOPE-DOC 100



Product introduction

- Telemedicine, and national health examination
- Portable, easy operation with high resolution
- Multiple functions, photo/video model
- USB interface for easy data communication

Dino Digital Video Otoscope

Product introduction

Connects to computer via USB
Allows the users to view and save high resolution images of Tympanic membrane and ear canal

- Save pathological images
- Super video record, able to set timer recording
- Hearing aid dispenser can do online remote diagnosis with the playback function.
- Magnify the ear canal with 200 times. The patients can also see their canals clearly.
- Handle focus

Brand: Dion-Lite
Accessory: Driver CD×1pc, Digital Video, Otoscopex1pc, Pouch×1pc



Product type



0302041-001
Model: AM311HT
Resolution: 300,000 pixel
Light source: LED ×4



0302038
Model: AM4113-EHT
Resolution: 1,300,000 pixel
Light source: LED ×6



Digital video otoscope magnifying Lens

Accessories of digital video otoscope

Model : 4# -0302014
3# -0302015

Product details



C132001-001
C132001-002
C132001-003

Dosing spoon

The dosing spoons enable an appropriate proper impression material.

Color: C132001-001 White
C132001-002 Green
C132001-003 Yellow



0301009

Mixing slab

Used with spatula for mixing ear impression

Dimension: 205x150mm



0301010

Mixing spatula

Used with slab for mixing ear impression material



0302024

Head light

For ENT doctors and dispensers checking ear canal

LED light: 1W Power supply: AAA×2

Battery not included



0902026

GSS-7 Tuning fork

Aluminum alloy tuning forks with single frequency pure tone

Brand: Soundlink
Type: 8192Hz ×1 4096Hz ×1 2048Hz ×1
1024Hz ×1 512Hz ×1 256Hz ×1
128Hz ×1



0303001

Listening tube

To test BTE,ITE,ITC,CIC

Material: Plastic Tube length: 390mm
Components: stetoclip,adaptor,and two eartips



0303002
0303002-005

Stetoclip with couple

To test BTE,ITE,ITC,CIC

Material: Plastic
Tube length: 390mm
Package: 25pcs/carton 0303002-005
Single package 0303002



0303004

Stetoclip with damper

To test BTE,ITE,ITC,CIC. With damper

Material: Plastic
Tube length: 390mm Plastic: 25pcs/carton



0303003

Adapter for listening tube

Accessory of stetoclip

Material: Silicone Size: Φ10.5x20mm



0701053-001
0701054-001
0701055-001

Insert foam tip

With smooth and comfortable material,probed for fit testing.
Certificate:RoHS, REACH, biological compatibility

Material: PU
Shape: Cylinder
Package:50pcs/bag
Use:Disposable
Type: 0701053-001 M Color: Yellow Length: 13mm Dia: 13mm
0701054-001 S Color: Beige Length: 13mm Dia: 10mm
0701055-001 L Color: Yellow Length: 13mm Dia: 18mm



0304003

Impression foam - XL

Perfectly conform to ear canals for accurate placement and protection. With thread

Colour: Blue
Dimension: Φ13x6mm Package: 50pcs/pack



0304004

Impression foam - L

Perfectly conform to ear canals for accurate placement and protection. With thread

Colour: Red
Dimension: Φ10.5x9mm Package: 50pcs/pack



0304005

Impression foam - M

Perfectly conform to ear canals for accurate placement and protection. With thread

Colour: Grey
Dimension: Φ9x7mm Package: 50pcs/pack



0304006

Impression foam - S

Perfectly conform to ear canals for accurate placement and protection. With thread

Colour: Yellow
Dimension: Φ7x6mm Package: 50pcs/pack



0505005

Tweezers bayonet shape

Used for picking up impression foam

Material: Metal Length: 160mm



0505011

Stainless steel round tweezers

Used for picking up impression foam round tip design to ensure safety

Length: 140mm



0301015

Scissors-straight pointed

Professional medical scissors with sharp shear edge for cutting sound tubes, elbows, etc.

Material: Stainless steel

Net weight: 60g

Length: 156mm



0301016-001

Tubing expander(Straight type-long)

Stainless steel tubing expander pliers with improved two-way lever spring

Model: 549-S

Type: Straight type-long



0301016-004

Tubing expander(Straight type-short)

Stainless steel tubing expander pliers with improved two-way lever spring

Model: 549-S

Type: Straight type-short



0301016-002

Tubing expander(45°Angeld)

Stainless steel tubing expander pliers with improved two-way lever spring

Model: 549-A

Type: 45°Angeld



0505012-001

Micro Ear Forceps(Alligator 8mm)

Delicate scissors action stainless steel forceps allow for greater precision to extract cerumen from the ear canal.

Dimension: Alligator 8mm



0505012-002

Micro Ear Forceps(Alligator 4.5mm)

Delicate scissors action stainless steel forceps allow for greater precision to extract cerumen from the ear canal.

Dimension: Alligator 4.5mm



0505012-003

Micro Ear Forceps(Cupped 4.5mm)

Delicate scissors action stainless steel forceps allow for greater precision to extract cerumen from the ear canal.

Dimension: Cupped 4.5mm



0902016

Sound Level Meter (1305A)

Features a range from 35 to 130dB, A&C frequency weightings, fast&slow dynamic characteristic models AC and DC signal output and back-screw for tripod fixture

Measurement Range: 35~130dB

Frequency response: 31.5Hz~8KHz



0902016-001

Sound Level Meter (Noise Measurement Tester)

Sound intensity testing tool

Measurement Range: 35~130dB

Frequency response: 31.5Hz~8KHz

MARKETING & DISPLAY

Focus on hearing aid accessories
Professional market display for your
hearing aids

Silicone ear model

Made of silicone
Ideal for shop window display and a practice tool for ear
impression taking
Other color and size is available on request
Either left or right is available

Dimension: 65x45x35 (H)mm



Blue

0101069-011 (R)
0101069-005 (L)



Green

0101069-009 (R)
0101069-003 (L)



Beige

0101069-012 (R)
0101069-006 (L)



Red

0101069-010 (R)
0101069-004 (L)



Transparent

0101069-014 (R)
0101069-008 (L)



Yellow

0101069-013 (R)
0101069-007 (L)



White

0101069-016 (R)
0101069-002 (L)



0103003-001
0103003-003

Cylindrical display (Single ear)

Provide unique and attractive visual effect for hearing aids!

Material : PC+ABS
Dimension : φ60x130 mm
Color : 0103003-001 Black
0103003-003 White

Color of ear model can be changed



0101033

Ear display (Single ear)

Provide unique and attractive visual effect to demonstrate the BTE, ITE, ITC and CIC

Dimension: 130x90x15 (H)mm
Content : One silicone ear model

Color of ear model can be changed



0101014

Ear display (Head style)

A display shelf in similar size of human head
Provide unique and attractive visual effect to demonstrate the BTE, ITE, ITC and CIC

Dimension: 190x150x15 (H)mm
Content : One silicone ear model

Color of ear model can be changed



0101014-012

Thicken Ear display (Frosted Head style)

Provide unique and attractive visual effect for hearing aids

Material: Acrylic
Dimension: 155x195x50(H)mm
Content: One pair silicone ear models

Color of ear model can be changed



0101014-011

Thicken Ear display (Clear Head style)

Provide unique and attractive visual effect for hearing aids

Material: Acrylic
Dimension: 155x195x50(H)mm
Content: One pair silicone ear models

Color of ear model can be changed



0101019-013

Ear display (Four ears)

Provide unique and attractive visual effect for hearing aids

Material: Acrylic
Dimension: 250x110x8 (H)mm
Display Color : Black
Content : Four silicone ear models

Color of ear model can be changed



0101019-011

Ear display (Four ears)

Provide unique and attractive visual effect for hearing aids

Material: Acrylic
Dimension: 250x110x8 (H)mm
Display Color : Transparent
Content : Four silicone ear models

Color of ear model can be changed



0101019

Ear display (Four ears)

Provide unique and attractive visual effect to demonstrate the BTE, ITE, ITC and CIC

Material: Acrylic
Dimension: 180x120x15 (H)mm
Content : Four silicone ear models

Color of ear model can be changed



0101031

Ear display (Five ears)

Provide unique and attractive visual effect to demonstrate the BTE, ITE, ITC and CIC

Material: Acrylic
Dimension: 250x180x15 (H)mm
Content : Five silicone ear models

Color of ear model can be changed



0101064

Ear display (Six ears)

Provide unique and attractive visual effect to demonstrate the BTE, ITE, ITC and CIC

Material: Acrylic
Dimension: 312x100x15 (H)mm
Content : Six silicone ear models

Color of ear model can be changed



0101065
0101065-001

Acrylic ear display

Imitated real ear canal of human beings.
Ideal to show how to wear earmolds and hearing instruments, and the position of the auditory canal. Suitable for practicing ear impression-taking

Material: Clear acrylic
Dimension: 62x42x45 (H)mm
Type: right ear-0101065 left ear -0101065-001

Hearing aid display

Provide unique and attractive visual effect to demonstrate BTE, ITE, ITC and CIC

Material: Acrylic
Dimension: 33x33x30 (H)mm

Colorful display

LED lights illuminate, rotate 360 degrees
Provide unique and attractive visual effect to demonstrate BTE, ITE, ITC and CIC

Power supply: AAx3 or DC 4.5V, 0.35A
Dimension: Φ95x40mm

Battery not included

Earmold display

High clear PE elastic film. Perfect to display earmold inside

Dimension: 0101073-001-Black, Out size: 108*128mm, Inner size: 92*92mm
0101073-002-White, Out size: 118*118mm, Inner size: 98mm
0101073-003-White, Out size: φ128mm, Inner size: φ104mm
A060301-128-White, Out size: 140*140mm, Inner size: 117*117mm
A060301-134-Black, Out size: 140*140mm, Inner size: 117*117mm

Crystal inner Carving Ear Anatomy Diagram

Acrylic 3D stereoscopic inner carving ear anatomy diagram allows the patient to clearly understand the structure of the human ear, it is convenient for the audiologist to analyze the condition of the ear for patients. This product can also be used as a window display.

Dimension: 198x50x98 mm- L
148x50x78 mm- S

OEM size and pattern are available on request



0102003

Ear anatomy model-L

Anatomical display model demonstrate outer, middle and inner ear

Material: Plastic
Dimension: 450x180x180(H)mm



0102012

Ear anatomy model-M

Anatomical display model demonstrate outer, middle and inner ear

Material: Plastic
Dimension: 290x130x180 (H)mm



0102010

Ear anatomy model-S

Anatomical display model demonstrate outer, middle and inner ear

Material: Plastic
Dimension: 115x75x105 (H)mm



0102005

Inner ear anatomy model

A professional model to demonstrate the organ of human inner ear, consisting of semicircular canal vestibule, cochlea and acoustic nerve

Material: Plastic
Dimension: 250x120x250 (H)mm



0102011-001

Ear anatomy diagram

Wall-mounted diagram to illustrate human ear anatomy

Dimension: 820x530 mm
Language: English



0104004-001

0104004-002

Hearing aid case- I

Small case for storing hearing aids

Material: Liner foam+ABS
Color: black / white
Dimension: 60x63x30 (H)mm
Package: 100pcs/carton



0104005-001

0104005-002

Hearing aid case

Small case for storing hearing aids

Color: 0104005-001 Black
0104005-002 White
Dimension: 70x75x27mm



0104005-003

0104005-004

Hearing aid case

Small case for storing hearing aids

Color: 0104005-003 Black
0104005-004 White
Dimension: 85x65x28mm



0104005-005

0104005-006

Hearing aid case

Small case for storing hearing aids

Color: 0104005-005 Black
0104005-006 White
Dimension: 65x65x30mm



0104005-007

0104005-012

Hearing aid case

Small case for storing hearing aids

Color: 0104005-007 White
0104005-012 Black
Dimension: 65x75x30mm



0104030-014

Hearing aid case

Small case for storing hearing aids

Material: Aluminum alloy
Color: Black
Dimension: φ80x35mm



0104014

Hearing aid case- II

Small case for storing hearing aids

Material: ABS
Dimension: 92x70x30 (H)mm
Package: 60pcs/carton



0104028-001

Hearing aid case

Small case for storing customized hearing aids

Material: ABS
Dimension: 60x77x29 (H)mm



0104029-001

Hearing aid case

Small case for storing hearing aids

Material: ABS
Dimension: 80x65x28 (H)mm



0104002

Presentation case- I

A delicate case to display hearing aids and accessories such as BTE, earmolds and customized hearing aids

Material: Aluminum alloy surface and plastic inside liner
Dimension: 255x170x40 (H)mm



0104003-005

Presentation case- II

A delicate case to display hearing aids and accessories such as BTE, earmolds and customized hearing aids
Without any samples inside

Material: aluminum alloy surface
Lining: upper cover- sponge, bottom surface-EVA
Dimension: 335x245x40 (H)mm



0104006

Hearing aid pouch

Handy storage for hearing aids

Material: PU
Dimension: 87x80mm



0104031

Hearing aid case

Small case for storing IEM, hearing aids, earmolds or hearing protectors

Material: EVA
Dimension: φ80x35 mm



0104018

Hearing aid case

A small case for storing hearing aids, earmolds and hearing protectors

Material: Mutispandex
Dimension: 96x70x33 (H)mm



0104022

Earphone case

To store hearing aids, earmolds and hearing protectors

Material: Oxford
Dimension: 115x85x50 (H)mm



0604011-002

Double-deck compartmental storage box

Easily taking out and holding hearing aid accessories

Material: PS
Dimension: 183x130x52 (H)mm



0604011-001

Nine drawer-style storage box

Easily taking out and holding hearing aid accessories

Material: PS
Dimension: 162x100x78 (H)mm

- 0903002-001
- 0903002-002
- 0903002-004

**Digital battery tester**

Small and easy-to-handle battery tester with LCD will quickly indicate the current power level with an ascending bar. Applicable sizes(675,13,312,10) with battery drawer and chain for attaching the key ring

Material: ABS
Color: black/red/blue
Package: 100pcs/carton

Certificate: CE, FCC



0903003

Activair battery tester

A large, easy-to-read indicator shows whether battery is good or bad, separate size slots for each type of batteries.

Material: ABS
Color: black
Package: 65pcs/carton

0904024

**Twin pack battery case**

A convenient battery case with snap-cover accommodating two 675 type cells or smaller
Non-breakable plastic construction with special loop for attaching to key chain or battery tester
Ideal for a giveaway with new hearing aid

Material: PP
Package: 50pcs/bag
Dimension: 36x18x7mm

EARMOLD LAB MATERIAL & EQUIPMENT

Supplying complete product
for earmold lab

NICE-FIT Series



500g 250g

Nice-fit hard ear mold material

For hard earmold products with UV light polymerisation
Brand: Nice-fit

Capacity: Transparent 500g 0402011-002
 Transparent 250g 0402012




0402008-002

Nice-fit shell material

Imported from Germany, used to make custom hearing aid and IEM shells.
Liquid UV-curing materials, high finished strength

Capacity: Transparent 600g




0403018
0403007
0403008



Nice-fit fast solidified silicone material

Perfect choice for making hearing aid soft earmolds
With good flexibility, high tensile strength, and comfortable wearing feeling
Solidified at room temperature

Brand: Nice-fit
Package: 2x25ml/pc, 8pcs/box
Base: square
Type: 25 shore rose 0403018
40 shore clear 0403007
40 shore rose 0403008



0403012
0403019
0403020



Nice-fit slow solidified silicone material

Perfect choice for making hearing aid soft earmolds
With good flexibility, high tensile strength, and comfortable wearing feeling
Solidified at room temperature

Brand: Nice-fit
Package: 2x25ml/pc, 8pcs/box
Base: square
Type: 25 shore rose 0403012
40 shore clear 0403019
40 shore rose 0403020



0403021

Nice-fit colorful slow solidified silicone material

Perfect choice for making hearing aid soft earmolds
With good flexibility, high tensile strength, and comfortable wearing feeling
Solidified at room temperature

Brand: Nice-fit
Package: 2x25ml/pc, 8pcs/box
Color: red, yellow, blue, green, orange, black, white, beige
Base: square
Type : 40 shore



Nice-fit Agar

Used for earmold duplicating

Brand: Nice-fit
Package: 0406005 1kg
0406004-004 6kg

Impression wax

Used to refine ear impression before making earmold

Brand: Nice-fit
Package: 500g/pack
Color: Rose



0408007

Polishing Paste(Rough)

A product for polishing hearing aid earmold and custom IEM shell

Brand: Nice-fit
Weight: 300g



0407002

Polishing Paste(Smooth)

Made In Germany, necessary material for hearing aid earmold and custom IEM shell polishing

Brand: Dreve
Weight: 900g



0407001-001



0405008-001

LP - Cleaner

Used for cleaning earmolds and ITE/ITC/CIC shells

Brand: Nice-fit
Package: 1L



Dreve shell material

Hard, thin shell-forming material(0.5-1.2mm can be regulated by polymerisation time). Identical wall thickness, especially for IEM hearing systems.

Brand: Dreve
Capacity: Transparent 500g-0404022-001
Transparent red 500g-0404022-002
Transparent blue 500g-0404022-003
Beige 500g-0404022-004
Transparent red 100g-0404022-005
Transparent blue 100g-0404022-006
Transparent dark 500g-0404022-007



0402011-001
0402012-001
0402013-001

Hard earmold material

Used for making hard earmolds , IEM shell , etc.

Brand: Dreve
Package: 0402011-001-500g- Transparent
0402012-001-250g- Transparent
0402013-001-500g- Light red

0403023-011
0403023-029




Silicon material Biopor AB

Used for making hearing aid soft earmolds
Can be used with manual inject pistol (0301004)

Brand : Dreve
Hardness: 25 shore
Specification: 50mlx8,
Item: 0403023-011 - Reddish transparent
0403023-029 - Transparent

Silicon material Biopor AB

0403023-001
0403023-021
0403023-027





Used for making soft earmold, sleeping plugs and swimming plugs

Brand: Dreve
Hardness : 40 shore
Specification: 50mlx8, round
Item: 0403023-001 - Reddish Transparent
0403023-021 - Transparent
0403023-027 - Blue Transparent

Hardness : 60shore
Item: 0403023-020 - Reddish transparent

Silicone material Biopor AB light

0403023-007
0403023-008





Used for making super soft earmold , sleeping plugs and swimming plugs

Brand: Dreve
Hardness: 25 shore
Specification: 50mlx8, round
Item: 0403023-007 - Bright green
0403023-008 - Purple

Silicon material Biopore AB Xtreme



0403023-022



Used for making soft earmold, sleeping plugs and swimming plugs



Brand: Dreve
Hardness: 25 shore
Specification: 50mlx8, round
Item: 0403023-022 - Transparent

Silicon material Biopore AB Xtreme



Used for making super soft earmold , sleeping plugs and swimming plugs

Brand: Dreve
Hardness: 40 shore
Specification: 50mlx8, round
Item: 0403023-017 - Transparent brown
0403023-018 - Transparent
0403023-019 - Reddish transparent
0403023-023 - Black
0403023-025 - White
0403023-026 - Grey



Silicon material Biopore AB Xtreme


Used for making super soft earmold , sleeping plugs and swimming plugs

Brand: Dreve
Hardness: 60 shore
Specification: 50mlx8, round
Item: 0403023-012 - Black
0403023-013 - Red
0403023-014 - Blue
0403023-015 - White
0403023-016 - Yellow

Swimming silicone material

Made in Germany, used for making customized swimming plugs

Brand: Detax
Specification: 50ml*8、 square
Hardness: 35shore
Color: Red x4 Blue x4




0403013



Direct molding silicone material

Made in Germany, used for making Soft earmold of BTE hearing aid

Brand: Detax
Specification: 50ml*8、 square
Hardness: 60shore
Color: Transparent red




0403014



Colorful swimming silicone material

Made in Germany, used for making Soft earmold of BTE hearing aid

Brand: Detax
Specification: 50ml*8、 square
Hardness: 35shore
Color: Red x2 Blue x2
Orange x2 Yellow x1
Green x1



0403022-001



0503002-001
Single earmold lab bench
HIGH EFFICIENCY LOW NOISE



Specification

- Operating desk
- Adjustable lighting
- Shelf
- Pneumatic gun
- Marble
- Dust drawer
- Subplate
- Suction nozzle
- Board to support hands
- Standby socket
- Silent vacuum system

Dimension


- L : 120 cm
- W: 60 cm
- H : 80 cm

Available colors



(The pictures are only for reference)

0503002-002
Double Earmold lab bench
HIGH EFFICIENCY LOW NOISE




Vacuum cleaner

- Absorbed amount 18-27 L/S adjustable
- Working noise < 65dB

Dimension

- L : 197 cm
- W: 60 cm
- H : 85cm

Available colors



(The pictures are only for reference)



0103007

Metal holder with clamp

For fixing the earmold during varnishing

Net weight: 60g
Material: Copper



0501017

Turning motor

Used for varnishing purpose
The rotation ensures that the lacquer evenly during the curing process

Power supply: AAA x 2
Dimension: 80x58x37(H) mm
Packing details: Single box + carton



0501016-005

LED curing unit

Used for curing 1-2 pcs of the earmold or shell of hearing aid & IEM one time. The LED light source is extremely durable and ensures a gentle polymerization of UV-material.

Brand: Dreve
Weight: 1.5 kg
Dimension: 205 x 205 x 255(D)mm
Voltage: 100-240V/50-60Hz



0501023-001
0501023-002

UV-LED curing unit

With LCD display, controlled by MCU technology.
High strength lightcuring time can be set customized.
Quickly setting red wax, lacquer and resin material for producing and maintain hearing aids& earmolds.

Brand: Soundlink
Model: GH05-I 0501023-001
Voltage: DC 24V/4A
Power: ≤60W
Weight: 3.5kg

Model: GH05-II 0501023-002
Voltage: DC 36V/4A
Power: ≤160W
Weight: 7.5kg



0503027-009

3D Nitrogen UV-curing unit

The machine is a light polymerization chamber for the hardening/curing resins. Two photoflash lamps at the bottom of the curing chamber operate in a work mode with a frequency of 10 flashes per second.

Voltage: AC220V/50Hz
Flashed lamp power: 100wx2
Curing chamber: 120x120x50(H) mm
Outer dimension: 310x310x140(H) mm
Capacity: 8pcs ITC/6pcs ITE
Spectral Distribution: 280-700nm, Max 400-500nm
Flash Rate: 10 flashes per second
Digital Time: Adjustable from 1up to 9999 flashes



0501009

Point lights UV-curing unit

A hand unit with light conductor for repairing hearing aid shell and making custom hearing aid vent

Weight: 2.2kg
Power: 75W
Voltage: 110-220V
Dimension: 200x140x135(H)mm



0501016-001

Open-ended UV-curing unit

Made in Germany
The unit has got a light source and mirrored interior surfaces with assure a safe polymerisation. With the battery operated turning motor, you obtain an even and completely cured surface.

Brand: Dreve
Weight: 1.7kg
Light source: 9Wx 2
Dimension: 240x95x100(H)mm
Light wave spectrum: 315-400nm
Power consumption: 35W(230V/50Hz)



0501009-002

Point lights UV-curing unit

Hand-held design, multi-time control options, efficient curing, especially suitable for partial varnish

UV light lamp



0501007
0501010
0501005
0501006
0501015

Replacement lamps for UV unit

Power : 0501007 9W 2Pin
0501026 9W 2Pin Osram
0501010 9W 2Pin Philips
0501005 11W 2Pin Philips
0501006 18W 4Pin Philips
0501015 36W 4Pin Philips

Point ling halogen lamp



0501014

Fit with light UV-curing unit

Power: 75W
Brand: Philips
Voltage: 12V

Hand drill



0502011-001

Brushless, high power, high torque unit
Features selection for forward and reverse direction
Safety overload protection, Torque feedback circuitry
With foot control model

Rotate speed: 0-50000r/min
Power: 230W
Voltage: 220V
Rated current: 2A
Max torque: 7.8N.cm

Hand drill



0502011-002

With brush hand drill, high power, high torque unit
Features selection for forward and reverse direction
Safety overload protection, Torque feedback circuitry
With foot control model

Rotate speed: 0-35000r/min
Power: 230W
Voltage: 220V
Rated current: 2A
Max.torque: 3.0Ncm



0503005-001

Electric heating knife

Used with the repairing wax to quickly and effectively repair the ear impressions

Power Consumption: 20W
Temperature Range: 50-200 °C
Power: AC100~240V
Dimension: 120x150x70mm



0503005-002

Electric heating knife

Originally imported from Germany, used to repair defective ear impressions. The temperature of the molten wax can be adjusted, and the heating state is displayed by the LED indicator

Heating power consumption: 7.5W
Machine power consumption: 10W
Temperature Range: 50-350 °C
Power: 220V/50HZ
Dimension: 100x65x55mm



0503034-001

AB glue machine

Used for mixing AB glue material

Brand: Dreve
Weight: 10kg
Voltage: 220V/50Hz
Power: 280W
Dimension: 405*225*319mm



0503034-002

AB glue charging barrel

Used for AB glue material

Brand: Dreve
volume: 1.7Lx2



0503034-003

Cable

Used for the connection between the AB glue machine and the charging barrel

Brand: Dreve



0503022

Duplicating Device

With the computer program and the unique heating pattern of patent technology, to automatically liquefy and maintain reversible duplicating gel at the correct temperature for duplication

Voltage: 220V/50Hz Capacity: 8kg
N.W.: 32kg Motor speed: 40rpm
Heating power: 800W×2 Dimension: 440x360x610(H) mm

Duplicating Device

High quality machine, LCD display
With thermostatic control steady functions, easy maintain.
Perfect machine for duplicating device

Brand: Simed Weight: 21kg
Heating Power: 400W Total Power: 600W
Power supply: 220V/50Hz Dimension: 350x370x430(H) mm

Wax melter

For wax melting. Easy to use.

Power: 500W
Weight: 0.92 kg
Power: 220V/50Hz
Dimension: 165x170x85mm

Wax melter

The device is suitable for temperature of approx. 10 °C (50 °F) to 100 °C (212 °F). The temperature control is infinitely variable. The container holds approx. 500 ml paraffin wax.

Brand: Dreve
Weight: 800g
Power: 100W/220-230V/50Hz
Dimension: 200x176x176mm

Digital Thermometer

Used to test temperature of wax or agar.
Good partner with wax melter(0506004)

Temperature range: -150 C~1300 C



0503006

Pressure pot

Made in Germany
Pressure polymerisation pot for silicone earmolds in water, no compressed air bubble and electricity supply required



0501012

Pressure pot

Made in Germany. Pressure polymerization unit with range up to 95 °C for one brass flask including clamp or four large rubber embedding molds.

Brand: Dreve
Model: Polymax-1
Weight: 10kg
Power:450W
Voltage: 220V/50Hz
Pressure range:3-5bar
Pressure pot: 160*110mm

Polishing and Grinding Machine- I

Perfect for modifying earmolds and hearing aids shells.

Brand: Soundlink
Weight: 36kg
Voltage: 220V/50Hz
Rotate speed: 3000r/min
Dimension: 600x500x460(H) mm



0503013-001

Grinding motor

Used to grinding hollow ITE/ITC/CIC shell or hard earmold

Brand: Soundlink
Voltage: 220V/50Hz
Power consumption : 100
Rotation rate: 2920r/min
Dimension: 160x140x153(H) mm



0501011



0506002



0506002-001



0301011



0504016

Plastic alcohol bottle

Used for serving alcohol
Liquid from the bottle can be sucked out by hand-operated pumps.

Cubage: 200ml



0504007

Cover plates

For the production of ITE/ITC/CIC shells, in order to prevent polymerization from above

Dimension: 44mm

Material: Acrylic



0504005-001

Antistatic wrist strap

Antistatic for making and fixing hearing aids with rope or without rope
Adjustable at will, rinsing resisting

IMQ: ±5%

Current limiting resistor

Electro-scattering time 0.1sec

Safe and reliable



0504003

Cutting knife

Used for cutting impression models



0504010



0504010-001

Blade

Used for cutting fine places of earmolds and hearing aids

Type: 10pcs /bag



0504009

Insulation paper

Used for insulation between the hearing aid components

Type: 100girds/paper



0703050-002

Tubing pullers

Used for pulling the tubing of soft earmold



Antistatic plastic zipper bags

Antistatic plastic bags with zipper for storing small parts and hearing aids

Dimension: 170×110mm - 0505010

170×240mm - 0505010-003

70×70mm - 0505010-004



0505001

Antistatic pointed tweezers

Used with anti-electrostatic wrist

Used for picking up small hearing accessories

Length: 120mm

Material: Plastic

Item: 0505001-pointed

0505002-flat

0505003-small flat



0505002



0505003



0505004

Stainless steel pointed tweezers

Used for picking up small hearing accessories

Length: 115mm

Material: Metal



0505005

Tweezers bayonet shape

Used for picking up small hearing accessories

Length: 160mm

Material: Metal



0505006

Tweezers flat shape

Used for picking up small hearing accessories

Length: 120mm

Material: Metal



0505011

Stainless steel round tweezers

Used for picking up small hearing accessories, round tip design to ensure safety



0504017

Glass bar

Used for wiping soft earmold material on the bottom of gel



0601005-002

Stereoscopic microscope PXS1020

Used for welding the hearing aid components

Dimension: 380×290×230 (H)mm
Distance between both eyes: 55-75mm
Weight: 4.2kg
Magnify: 10X,20X
Working Distance: 65mm
The workbench Diameter: 75mm



0602008-003

Weller Welding Unit WT1014

Used for manufacturing and fixing hearing aids, hyperthermal LCD display, anti-static function

Power: 90W
Temperature: 50°C-550°C
Model: WT1014
Voltage: 220V/50HZ



0602002 0602003

Welding head

The head of welding unit, connect with weller welding

Model: LT1S -0602002
LT1 -0602003
ETS -0602004
LTKN -0602006



0602012

Cleaning ball

Used with weller welding WD 1000 (0602010)



0601002

Bench magnifier

Hands-free magnified lighting for close, precision work
Attached with clamp
High quality desktop design
Portable pedestal base for hearing aid and earmold maintaining

Diameter of lens: 127mm
Color: white
Weight: 3.2kg
Magnify: 5X
Voltage: 200V



0604020
0604020-001

Lead-free solder wire

Weight: 250g/roll
Diameter: 0.4mm-0604020
0.3mm-0604020-001




0601003

Magnifier lamp

Optimal for control and soldering works.
The magnifier lamp disposes of a 3 dioptic optician lens
Attached with clamp

Magnify: 5X
Votage: 200V
Weight: 3.2kg
Diameter of lens: 127mm



Silver plated litz wire

Made from silver plated copper wire with two different layers of enamel
Used for connecting receiver, microphone and other accessories of hearing aids

Length: 10m
Specification: 10m/bag
Color: Red, Green, Beige, Black
Model: Red-0704005
Green-0704007
Beige-0704023-002
Yellow-0704006
Blue-0704023-001
Black-0704023-003



Soldering band

0604006 Used for welding hearing aids and IEM



Cutting plier

0604008 Side cutting nipper for miniature work

Length: 125mm Weight: 60g



Blue tack

0604002-001 For hearing aids repairing, ensure the surface clean and dry

Weight: 75g



Diagonal cutting plier

0604021 Used for cutting all kinds of wires



Casting ring

0504001 Used for making ear model with removable bottom part
0504015 Reusable forms in two different sizes
Transparent allows UV light to penetrate

Dimension: Φ 50x40mm- 0504001 Φ 41x37mm- 0504015
Material: Acrylic / PC Color: transparent



Impression lifter

0504002 Produced by stainless steel anti-skid design, used for catching impression diip with wax

Length: 175mm



Solder wire

0604007 For hearing aid welding



Waxdip filter

0506005 Used for ear impression waxdip, easy operation, suitable for learner to handle



Cutting-wire tweezers

0604001-001 Used for cutting the wires of electrical components



Wild-mouth bottle

0405020-002 Dark brown, used for storing resin material or lacquer

Volume: 100ml Dimension: Φ 70x63mm
Material: glass bottle, alumite cover

DRILLING & GRINDING



Ball burs

Popular assortment of steel, ball shaped burs for hard earmold sound outlet
Single cut

Aggressively cut & shape

- 0.5mm-1101022 2.3mm-1101026
- 1.0mm-1101023 2.4mm-1101027
- 1.2mm-1101001 2.5mm-1101028
- 1.3mm-1101024 2.7mm-1101005
- 1.6mm-1101002 3.0mm-1101029
- 1.8mm-1101003 3.1mm-1101006
- 2.0mm-1101025 3.5mm-1101020
- 2.1mm-1101004




1101007

Oval cutter

Popular assortment of steel, oval shaped burs for hard earmold grinding.
Single cut

Aggressively cut & shape

5.0



1101008
1101009
1101017

Cutters conically-shaped

Popular assortment of steel, conically shaped burs for hard earmold grinding.
Single cut

Aggressively cut & shape

- 1101008-4.0mm
- 1101009-5.0mm
- 1101017-6.0mm




1101019
1101012

Cutters pear-shaped

Popular assortment of steel, pear shaped burs
Single cut

Aggressively cut & shape

- 1101019-4.0mm
- 1101012-5.0mm




1101014
1101013

Cutters conically-shaped

Popular assortment of steel, ball shaped burs
Single cut

Aggressively cut & shape

- 1101014-4.5mm
- 1101013-5.0mm



1101021
1101011


- 1101021-5.0mm
- 1101011-6.0mm



1101018

Tungsten steel grinding

Popular assortment of steel, conically shaped burs
Single cut
Aggressively cut & shape



1102001
1102002
1102003
1102004

Soft earmold ball burs

Double cut burs remove material for soft earmold sound outlet drill
1102001-1.6mm
1102002-2.1mm
1102003-2.3mm
1102004-3.1mm



1102020
1102005

Soft earmold oval burs

Double cut burs remove material for soft earmold drill
1102020-5.0mm
1102005-7.0mm



1102006

Soft earmold cylindric cutter

Double cut burs remove material for soft earmold drill

5.0




1102007

Soft earmold pear cutter

Double cut burs remove material for soft earmold drill

6.0



001 002 003 004 005

Soft earmold cutter system

Carbonize drill for soft earmold resin grinding

1102028-001-8mm
1102028-002-5mm
1102028-003-8mm
1102028-004-8mm
1102028-005-8mm



1102008

5.0

Cylindrical grinding caps holder

Used with cylindrical grinding caps for impression and earmold



1102009
1102010
1102022

Cylindrical grinding caps

Used with cylindrical grinding caps holder for impression and earmold

Cylindrical grinding head 80 (1102009)
Cylindrical grinding head 150 (1102010)
Cylindrical grinding head 320 (1102022)

Package: 5pcs/pack



1102023

5.0

Conical grinding caps holder

Used with conical grinding caps for impression and earmold



1102024
1102025
1102026

Conical grinding caps

Used with conical grinding caps holder for impression and earmold

Cylindrical grinding head 80 (1102024)
Cylindrical grinding head 150 (1102025)
Cylindrical grinding head 320 (1102026)

Package: 5pcs/pack



Rubber polisher

Grinding the surface of impression, earmold, hearing aid shell, both elaborate polishing and grinding

- Rubber grinding head C-2 (1102012)
- Rubber grinding head C-3 (1102013)
- Rubber grinding head M-1 (1102014)
- Rubber grinding head M-2 (1102015)
- Rubber grinding head M-3 (1102016)
- Rubber grinding head F-2 (1102017)
- Rubber grinding head F-3 (1102018)



1103003

Cylindrical sandpaper holder

Used with sandpaper



1103018
1103025
1103017

Sandpaper

Used with sandpaper holder for impression, hard earmold and hearing aid shell grinding

- Dimension: 50m/roll
- Grain: 120 (1103018)
- Grain: 220 (1103025)
- Grain: 320 (1103017)



1102019

Rubber polisher 28

Grinding the small part of impression, earmold, hearing aid shell



1103024
1103033
1103030

Cylindrical grinding caps

For grinding earmolds and hearing aid shells
Accessory of grinding motor

- Diameter: 115mm
- Grain: 80 (1103024)
- Grain: 120 (1103033)
- Grain: 220 (1103030)



1102011
1102011-001
1102011-002
1102011-003

Soft earmold punch

Operating without hand drill, punching for soft earmold easily and fast
Necessary tool for soft earmold drilling

- Diameter of hole: 3.0mm(1102011)
- 3.5mm(1102011-001)
- 3.2mm(1102011-002)
- 3.8mm(1102011-003)



1103041-001
1103041-002

Sandpaper

For grinding earmolds and hearing aid shells
Accessory of grinding motor

- Dimension: 280x230mm
- Grain: 220 (1103041-001)
- Grain: 120 (1103041-002)




1102021-001
1102021-002

Soft earmold punch(special for children)

Operating without hand drill, punching for soft earmold easily and fast
Necessary tool for soft earmold drilling
Specially fit for patients with thinner canal and children

- Diameter of hole: 2.0mm(1102021-001)
- 2.5mm(1102021-002)



1103034
1103005
1103006
1103007
1103008

Twist drill

Great tools for hard earmold sound outlet drill

- Twist drill-0.5mm (1103034)
- Twist drill-0.6mm (1103005)
- Twist drill-0.8mm (1103006)
- Twist drill-1.0mm (1103007)
- Twist drill-1.2mm (1103008)
- Twist drill-1.4mm (1103009)
- Twist drill-1.6mm (1103010)
- Twist drill-1.8mm (1103038)
- Twist drill-2.0mm (1103004)
- Twist drill-2.3mm (1103021)
- Twist drill-2.35mm (1103022)
- Twist drill-2.85mm (1103011)



1103002

Conical sandpaper holder

Used with sandpaper



1103009
1103010
1103038
1103004
1103021
1103022
1103011




1103020
1103020-001

Speed Trimmer

Used for cutting faceplate of customized hearing aids

All-in-one (1103020)
Split type (1103020-001)



Finishing stone

Used for grinding impression mold

Color: blue
Cone 10mm (1103019)
Cone 10.5mm (1103031)
Cone 13mm (1103039)
Cone 13.5mm (1103036)
Cylindrical 13.5mm (1103037)



1103026
1103027

Emery ball burs

Used for drilling hard earmolds and hearing aids

Material: emery
2.3mm (1103026)
3.5mm (1103027)



1104001

Drill base

Used to put drills in order

Material: POM
Net weight: 300g
Diameter: 90mm



1102029

Soft earmold arc scoop

Efficiently cutting the duplicated gel to get the perfect arc shape of soft earmolds

Material: Metal



1103013
1103013-001

Wheel buffing-thick

Made of high quality muslin with combed edges for buffing hard earmolds and hearing aid shells

Material: Muslin
Diameter: 100mm
Hole Diameter: ϕ 10mm(1103013)
 ϕ 4mm (1103013-001)



1103029
1103029-001

Wheel buffing-thin

Thick buffing, made of high quality muslin with combed edges, for buffing hard earmolds and hearing aid shells

Material: Muslin
Diameter: 100mm
Hole Diameter: ϕ 10mm(1103029)
 ϕ 4mm (1103029-001)



1103014

Bufs hand shank

Fits with soft buffs (1103015), detachable and convenient



1103015

Soft buffs

Fits on buffs hand shank (1103014) for cleaning dust of earmold and hearing aid shell, and buffing prehiminary

Package: 10 pcs per pack



1103032

1103035

Polishing buff

Polish and grind the earmold and hearing aid
The tip of buff is non-removable.

1103032-Polishing buff
1103035-Wheel brush

HEARING AID ACCESSORIES



0701028 ■
0701029 ■



2pin receiver

Body aids accessory,used with 2pin cord

NO.	Item	Specifications
1	Impedance(Ω)	220 \pm 20% at 1,000Hz
2	D.C.Resistance(Ω)	90 \pm 10%
3	Sensitivity(dB)	116 \pm 3dB at 1,000Hz
4	Nominal Input	1mW
5	Maximum Input	300mW
6	Operating Temperature(C)	-20~70 C
7	Material	Cap Al Finish ABS Housing
8	Color	Black
9	Package	50pcs/pack

0701002 ■
0701004 ■



3pin receiver

Body aids accessory,used with 3pin cord

NO.	Item	Specifications
1	Impedance(Ω)	220 \pm 20% at 1,000Hz
2	D.C.Resistance(Ω)	45 \pm 10%
3	Sensitivity(dB)	116 \pm 3dB at 1,000Hz
4	Nominal Input	1mW
5	Maximum Input	300mW
6	Operating Temperature(C)	-20~70 C
7	Material	Cap Al Finish ABS Housing
8	Color	Black
9	Package	50pcs/pack

0701006-001



2pin bone conduction connector

Connects to 2pin cord for Bone conduction headphone

Model: 2Pin
Dimension: 16.7 \times 24 \times 14 (H)mm

0701005-001



3pin bone conduction connector

Connects to 3pin cord for Bone conduction headphone

Model: 3Pin
Dimension: 16.7 \times 24 \times 14 (H)mm

Standard 2pin cord

Body aids accessory,used with 2pin receiver

0701006 ■

0701027 ■



NO.	Item	Specifications
1	Length	750mm
2	Intertwist Distance	5mm
3	Lay Direction	S
4	Needle	1.28/1.58
5	Spacing of Needle	3.8mm
6	Wire diameter	1.3-1.4mm
7	Material	PVC
8	Medium	Oxygen-free copper
9	Resistance(Ω)	0.6-0.8
10	Color	Beige
11	Package	50pcs/pack

High-quality 3pin cord

Body aids accessory,used with 3pin receiver

0701011 ■

0701023 ■



NO.	Item	Specifications
1	Length	750mm
2	Intertwist Distance	4mm
3	Lay Direction	S
4	Needle	1.18
5	Spacing of Needle	2mm
6	Wire diameter	1.3-1.35mm
7	Material	PVC
8	Medium	Copper wire
9	Resistance(Ω)	1.8-2.0
10	Color	Beige
11	Package	50pcs/pack

High-quality 2pin cord

Body aids accessory,used with 2pin receiver

0701008 ■

0701009 ■



NO.	Item	Specifications
1	Length	750mm
2	Intertwist Distance	5mm
3	Lay Direction	S
4	Needle	1.28/1.58
5	Spacing of Needle	3.8mm
6	Wire diameter	1.2-1.3mm
7	Material	PVC
8	Medium	Copper wire
9	Resistance(Ω)	0.8-1.0
10	Color	Beige
11	Package	50pcs/pack

2pin Y cord (normal)

Body aids accessory,used with 2pin receiver

0701012 ■

0701037 ■



NO.	Item	Specifications
1	Length	750mm
2	Intertwist Distance	5mm
3	Lay Direction	S
4	Needle	1.28/1.58
5	Spacing of Needle	3.8mm
6	Wire diameter	1.3-1.4mm
7	Material	PVC
8	Medium	Oxygen-free copper
9	Resistance(Ω)	0.8-1.0
10	Color	Beige
11	Package	25pcs/pack

Standard 3pin cord

Body aids accessory,used with 3pin receiver

0701010 ■

0701038 ■



NO.	Item	Specifications
1	Length	750mm
2	Intertwist Distance	4mm
3	Lay Direction	S
4	Needle	1.28/1.58
5	Spacing of Needle	2.0mm
6	Wire diameter	1.3-1.4mm
7	Material	PVC
8	Medium	Oxygen-free copper
9	Resistance(Ω)	0.6-0.8
10	Color	Beige
11	Package	50pcs/pack

2pin Y cord (advanced)

Body aids accessory,used with 2pin receiver

0701013 ■

0701014 ■



NO.	Item	Specifications
1	Length	750mm
2	Intertwist Distance	5mm
3	Lay Direction	S
4	Needle	1.28/1.58
5	Spacing of Needle	3.8mm
6	Wire diameter	1.2-1.3mm
7	Material	PVC
8	Medium	Copper wire
9	Resistance(Ω)	0.8-1.0
10	Color	Beige
11	Package	25pcs/pack



0701034 ■
0701035 ■

3pin Y cord

Body aids accessory,used with 3pin receiver

NO.	Item	Specifications
1	Length	750mm
2	Intertwist Distance	4mm
3	Lay Direction	S
4	Needle	1.18
5	Spacing of Needle	2.0mm
6	Wire diameter	1.3-1.4mm
7	Material	PVC
8	Medium	Oxygen-free copper
9	Resistance(Ω)	0.6-0.8
10	Color	Beige
11	Package	25pcs/pack



0701007

Single 2pin cord

Body aids accessory, used with 2pin receiver

Color: Beige
Length: 750mm
Package: 50pcs/pack



0701046-001
0701046-002

USB Cable

Body aids accessory,used with 2pin receiver

Color:white
Unilateral (0701046-001)
Bilateral (0701046-002)



0701050-001
0701050-002
0701050-003

Open fitting eartip

Replacement tips for hearing aid

Material: Silicone
Dimension: 6mm (0701050-001)
9mm (0701050-002)
12mm (0701050-003)



0701050-007
0701050-008
0701050-009

Open fitting eartip- II

Double-layer eartips to meet the needs of people with heavy hearing loss

Standard hardness: Large (0701050-007)
Medium (0701050-008)
Small (0701050-009)



0701050-010
0701050-011
0701050-012

Soft: Large (0701050-010)
Medium (0701050-011)
Small (0701050-012)



0701050-013
0701050-014
0701050-015

Open fitting eartip- III

Replacement for open fitting hearing aid.

Color: Transparent Black
Dimension: φ10mm (0701050-013)
φ8mm (0701050-014)
φ6mm (0701050-015)



0701050-017

Open fitting eartip- IV

Replacement for open fitting hearing aid.

Color: Transparent Black
Dimension: 14*11*9mm



0701050-018

Open fitting eartip- V

Replacement for open fitting hearing aid.

Color: Transparent
Dimension: φ8*7mm



Eartip- II

Replacement tips for Body Aids and BTE

Material: Silicone
 Package: 100pcs/pack
 Certificate: RoHS, Biocompatibility



Eartip with tubing- II

A combination of eartip,tubing connector and tubing (3.2×2.0mm) fit for BTE

Package: 100pcs/pack



PVC tubing with connector- II

A combination of eartip,tubing connector and tubing(3.2×2.0mm) fit for BTE

Package: 100pcs/pack



Earplug adaptor- II

For body-worn hearing aid receiver

Material: PP
 Package: 100 pcs/pack



Tubing connector- II

Standard hearing aid accessory for BTE.
 Used to connect eartip to tubing (3.2×2mm)

Material: PP
 Package: 100 pcs/pack



Silicone tubing

Elastic material easy to be stretched.
 Specially used for small ear canal of children

Material: Silicone
 Dimension: 3.0×1.8mm



Silicone tubing

Used to produce vent in customized hearing aids,
 can't be heated,when fixed by gule.

Material: Silicone
 Pack: 10m/bag
 Dimension: 1.4×1.0mm



PVC tubing

Standard hearing aid accessories for BTE and earmold

Material: PVC
 Package: 100m/roll
 Diameter: 3.0×2.0mm(0703001)
 3.5×2.0mm(0703008)
 3.3×2.0mm(0703053)
 3.2×2.0mm(0703028)
 1.5×2.5mm(0703046)
 1.0×2.0mm(0703003)
 4.0×2.5mm(0703022)



Moistureproof PVC tubing

Standard hearing aid accessories for BTE and earmold, with
 moistureproof design

Material: PVC
 Package: 100m/roll
 Diameter: 3.0×2.0mm



Moistureproof Bent tubing

Used for hard earmold

Specifications: 3.0×2.0mm(0703042)Moistureproof Bent tubing
 3.0×2.0mm(0703042-001)Moistureproof Prebent tubing



0703027
0703031
0703032
0703055-001

Bent tubing

Used for soft earmold, connects with earmold and BTE

Material: PVC

Package: 100pcs/pack

Specifications: 3.5×2.0mm(0703027)
3.2×2.0mm(0703032)
3.0×2.0mm(0703031)
2.5×1.5mm(0703055-001)



0703010

Earmold hook- III

Used for attaching tubing to earmolds

ID: 1.5mm

Material: Nylon

Package: 100pcs/pack



0703026
0703033
0703052
0703051-001
0703014-001

Prebent tubing

Used for soft earmold, connects with earmold and BTE

Material: PVC

Package: 100pcs/pack

Specifications: 3.5×2.0mm(0703026)
3.3×2.0mm(0703052)
3.2×2.0mm(0703033)
3.1×2.0mm(0703051-001)
4.0×2.0mm(0703014-001)



0703020

Tube lock- I

Used to match bent tubing (3.5×2.0mm) for soft earmold

Material: Nylon

Package: 100 pcs/pack



0703051
0703031-003

Bent tubing in oblique cut

Used for soft earmold, connects with earmold and BTE

Material: PVC

Specifications: 3.1×2.0mm(0703051)
3.0×2.0mm(0703031-003)



0703021-003



0703021-001

Tube lock- II

Used to match bent tubing (3.2×2.0mm) for soft earmold

Material: Copper

Package: 100pcs/pack



0703033-002
0703051-002
0703026-001

Prebent tubing in oblique cut

Used for soft earmold, connects with earmold and BTE

Material: PVC

Package: 100 pcs/pack

Specifications: 3.2×2.0mm(0703033-002)
3.1×2.0mm(0703051-002)
3.5×2.0mm(0703026-001)



0703021-002

Connecting Clip

Material: copper, gilded

Dia: 1mm



0703005
0703011

Earmold hook

Used for hard earmold

Material: PP

Package: 100pcs/pack



0703016



0703015



0703019

Receiver tubing

For producing and repairing hearing aids

ID: 1.3 mm(0703015) 1.1mm(0703016) 1.5mm(0703019)

Material: Viton

Package: 100 pcs/pack



0703017
0703029

Receiver tubing

For producing and repairing hearing aids

ID: 1.45 mm(0703017) 1.85 mm(0703029)
Material: Viton Package: 100 pcs/pack



0704026-003
0704026-002
0704026-005
0704026-004

HF4 Wax filter

Used for preventing ear wax effectively

Item No: HF4 Wax filter red with tool (16 pcs/pallet) - 0704026-003
HF4 Wax filter blue with tool (16 pcs/pallet)- 0704026-002
HF4 Wax filter red - 0704026-005
HF4 Wax filter blue- 0704026-004



0704001

Wax spring

Prevent earwax from entering into the tubing of ITE, ITC and CIC

Material: Copper Package: 100 pcs/pack




0704026-013 0704026-014

Wax filter

Used for preventing ear wax effectively

Model: 1.0mm White 0704026-013
1.2mm Grey 0704026-014



0704012
0704013

Wax guard

Used for customized hearing aids,connect with wax filter

Band: Resound Spec: 100pcs/bag
Model: Red 0704012 Blue 0704013



0704024-001

Wax filter (Nanocare)

Used for preventing ear wax effectively



0704014

Wax filter

Used for customized hearing aids,use with wax guard


Band: Resound Pack: 10pcs/plate 60pcs/box



0704026-001

Wax filter (Connexx)

Used for preventing ear wax effectively



0704015

Wax filter

Used for customized hearing aids,use with wax guard

Band: Resound Pack: 100pcs/bag



0704017 0704016

Cerumen protection series

Used for customized hearing aids, prevent earwax effectively

Red: Right ear Blue: Left ear
Pack:100pcs/bag
Model: Red cover clear base-0704016
Blue cover clear base-0704017
Red cover red base-0704018
Blue cover blue base-0704019



0704019

0704018

Venting Wire



Used to produce the vent in customized hearing aids

Brand: Dreve Pack: 10 pcs/bag
Model: 0.8mm-0704021-001 1.0mm-0704021-002
1.2mm-0704021-003 1.5mm-0704021-004

Pull out string




0704011-001

Used for taking out customized hearing aids

Total length: 43.5±0.1mm
Bead Diameter: 1.5±0.4mm String length: 35±5mm

Silver plated litz wire

Made from silver plated copper wire with two different layers of enamel
Used for connecting receiver, microphone and other accessories of hearing aids



Length: 10m
Specification: 10m/bag
Color: Red, Green, Beige, Black
Model: Red-0704005 Yellow-0704006
Green-0704007 Blue-0704023-001
Beige-0704023-002 Black-0704023-003

Length: 30mm
Specification: 500pcs/box

Dia of wire: 0.05mm*7
Color: Red, Green, Yellow, Blue
Model: Red-0704020-019
Green-0704020-021
Yellow-0704020-020



Dia of wire: 0.3mm
Color: Red, Green, Gold, Black
Model: Red-0704020-028
Green-0704020-029
Gold-0704020-030



0705



0706



0720

Knowles Receiver System

Standard hearing aid & IEM accessories from Knowles, several models available

CI	EF	EP	EH	FK	FH	ED		BK
Specification	Specification	Specification	Specification	Specification	Specification	Specification		Specification
CI-22955	EF-21937	EP-27110	EH-27232	FK-23451	FH-23281	ED-21913	ED-27962	BK-21600
CI-28267	EF-21942	EP-24074	EH-23149	FK-23466	FH-23853	ED-27236	ED-27926	BK-21606
CI-28352	EF-21958	EP-24075	EH-27822	FK-26260	FH-23851	ED-23146	ED-21932	BK-21613
CI-28329	EF-21943	EP-24123	EH-27827		FH-26664	ED-23162	ED-23147	BK-21610
CI-26697	EF-27505	EP-27104	EH-23043		FH-23373	ED-29457	ED-23619	BK-21968
CI-22764	EF-26356	EP-27107			FH-23375	ED-29689	ED-26805	BK-28507
		EP-27108			FH-26553	ED-28313	ED-23313	
		EP-24073			FH-23855	ED-23814	ED-27156	
					FH-23821	ED-30262	ED-27304	
					FH-23371	ED-27305	ED-27288	
					FH-23377	ED-27288	ED-27488	
						ED-21950		

Knowles Microphone

Standard hearing aid accessories from Knowles, several models available


EM		EK		EA	TM	FG
Specification		Specification		Specification	Specification	Specification
EM-25346	EM-25046	EK-223133	EK-25033	EA-221868	TM-23568	FG-23452
EM-23446	EM-24446		EK-23027	EA-21842	TM-23546	FG-23453
EM-24368	EM-23368	EK-23032	EK-23024	EA-23178	TM-24546	FG-23629
EM-23346	EM-23046	EK-25024			TM-26613	FG-23329
EM-26392	EM-23346	EK-23028				FG-26163
EM-26596	EM-25446	EK-26695				FG-23653
EM-24346	EM-23800	EK-23142				FG-60019

Knowles Damper

Standard hearing aid accessories from Knowles, several models available

Model Variations					
Color	Nom.Acoustic Resistance(Ω)	2.08mm Plug	1.78mm Screen	1.37mm Screen	1.12mm Screen
Grey	330	BF-1999	BF-3034	---	---
White	680	BF-1859	BF-1997	BF-7083	BF-1988
Brown	1000	BF-1860	BF-3035	BF-3292	BF-3320
Green	1500	BF-1861	BF-3038	BF-3039	BF-1991
Red	2200	BF-1921	BF-3037	BF-9778	BF-1995
Orange	3300	BF-1922	BF-3036	---	---
Yellow	4700	BF-1923		---	BF-3163

INDUSTRIAL NOISE PROTECTOR



CUSTOMIZED EARMOLDS

Currently more and more people pay more attention to the healthy of ear. Customized earmolds product, like IEM, Swimming earplug, Noise production, Sleeping earplug and phone earplug are welcomed in the market.

SLEEP EARPLUG



SWIMMING EARPLUGS



0403013

Swimming silicone material

Made in Germany, used for making customized swimming plugs

Brand: DETAX
Specification: 50ml*8、 square
Hardness: 35shore
Color: Red x4 Blue x4



0403022-001

Colorful swimming silicone material

Made in Germany, used for making Soft earmold of BTE hearing aid

Brand: DETAX
Specification: 50ml*8、 square
Hardness: 35shore
Color: Red x4 Blue x4
Orange x2 Yellow x1
Green x1



1102027

1102027-001

Soft earmold handle tool

Reusable used for making swimming earplug handle

Type: Arc 1102027
Cone 1102027-001



0205032
0205032-001
F2001004-003

Silicon rope

Waterproof silicone material, Floating in the water
Prevent from falling off

Type: 0205032 White dia 2mm
0205032-001 Black dia 2.2mm
F2001004-003 Blue dia 2.2mm



0205033
0205034

Cord lock

Used for R&L earplug in red and blue, used for 2mm white silicon rope

Material: ABS
Package: 10pcs/pack
Type: 0205033 Red
0205034 blue




F2001013-001

Rotating clamp

Can rotate 360°, fix the white dia 2mm silicon rope
connect swimming earplug and swimming cap

Material: ABS

Customized travel noise earplug





0901001

3M Ear plugs

Smooth,tapered surface.The soft foam material compresses then expands to conform to the shape of the ear canal.

Non-allergenic
Come by the pair

Brand: 3M NRR: 29dB
Material: PU Color: bright orange



0901011-003

3M Ear plugs

Smooth,tapered surface.The soft foam material compresses then expands to conform to the shape of the ear canal

Non-allergenic
Come by the pair with line,effectively prevent falling off

Brand: 3M NRR: 29db
Material: PU Color: bright green



0901002

3M Ear plugs with line

Smooth,tapered surface
The soft foam material compresses then expands to conform to the shape of the ear canal

Non-allergenic
Come by the pair with line,effectively prevent falling off

Brand: 3M NRR: 29dB
Material: PU Color: bright orange



0901004

3M noise reduction plugs with line

Ideal to protect hearing and reduce the noise
Come by the pair with line,effectively prevent falling off

Brand: 3M NRR: 23dB Material: Silicon



0901011-001

Honeywell earplug

Ideal to protect hearing and reduce the noise
Come by the pair

Brand: Honeywell Model: MAX-1



0901011-002

Honeywell earplug

Ideal to protect hearing and reduce the noise
Come by the pair with line, effectively prevent falling off

Brand: Honeywell Model: MAXLITE LPF-1



0901006

Silicone ear plugs

Ideal to protect hearing and reduce the noise

NRR: 25dB
Material: Silicone Size: φ20x5mm
2pcs/box



0901007

3M ear defenders

Used for protecting against noise
Constructed for comfort,soft foam filled cushions
Adjustable of the headband position(Head top,Head back and chin)

Brand: 3M NRR: 27dB Specification: 1427



0901008

3M ear defenders

Used for protecting against noise
Constructed for comfort,soft foam filled sponge

Brand: 3M NRR: 23dB Specification: 1425

IEM cable



0701047-019
0701047-020

This replacement cables fit all 2-Pin ES and AC Series Westone monitors as well as Universal Fit UM2 and UM3X models w/Removable Cables. It also fit other brands of musician monitors with two pin connectors.

Color: Silver
Pin connector diameter : 0.78 ± 0.01 mm
Plug cross section dimension: 4.6 mm * 2.8 mm (1.81 inch * 1.10 inch)



0701047-019
2PIN IEM Cable
Length: 127 cm (50 inch)



0701047-020
2PIN IEM Cable
Length: 162 cm (64 inch)



0701047-017
0701047-018



0701047-017
2PIN without MIC-Black



0701047-018
2PIN with MIC-Black



0701049-001
0701049-002



0701049-001
Memory wired
Length: 130 cm



0701049-002
Line curved
Length: 129 cm



0701049-010
0701049-011

0701049-010 MMCX IEM Cable -Black(Length: 166 cm)
0701049-011 MMCX IEM Cable-Silver(Length: 169 cm)

Socket

One of the IEM parts, used with the earphone cable

Package: 100 pcs/pack

Model: 0701041-001 H shaped 2Pin female socket
0701041-004 H shaped 2Pin female socket



Model: 0701048-001 2Pin Recessed female socket
0701048-002 MMCX female socket
0701048-003 2Pin female socket
0701048-004 2Pin female socket

Nice-fit hard earmold material

For hard earmold products with UV light polymerisation

Brand: Nice-fit

Capacity: Transparent	500g	0402011	<input type="checkbox"/>
Transparent	250g	0402012	<input type="checkbox"/>



0402011
0402012

Dreve IEM shell material

Hard, thin shell-forming material(0.5-1.2mm can be regulated by polymerisation time). Identical wall thickness, especially for IEM hearing systems.

Brand: Dreve


Capacity: Transparent	500g-0404022-001	<input type="checkbox"/>
Transparent red	500g-0404022-002	<input type="checkbox"/>
Transparent blue	500g-0404022-003	<input type="checkbox"/>
Beige	500g-0404022-004	<input type="checkbox"/>
Transparent red	100g-0404022-005	<input type="checkbox"/>
Transparent blue	100g-0404022-006	<input type="checkbox"/>
Transparent dark	500g-0404022-007	<input type="checkbox"/>



0404022-001-007

Silver plated litz wire

Made from silver plated copper wire with two different layers of enamel
Used for connecting receiver, microphone and other accessories of hearing aids



Length: 10m
Specification: 10m/bag
Color: Red, Green, Beige, Black
Model: Red-0704005 Yellow-0704006
 Green-0704007 Blue-0704023-001
 Beige-0704023-002 Black-0704023-003



Receiver for IEM



Complete series of IEM receiver, 100% imported from Knowles, quality assurance
Popular models: SR,RAB,ED,CI,BK,TWFK,DWFK,HODVTEC,WBTK, SWFK,GQ,HE,GK,GV

Length: 30mm
Specification: 500pcs/box



Dia of wire: 0.05mm*7
Color: Red, Green, Yellow, Blue
Model: Red-0704020-019
 Green-0704020-021
 Yellow-0704020-020



1105001

Damper Mounting Tool

Used for installing and dismounting damper of IEM

Length:36mm

Dia of wire: 0.3mm
Color: Red, Green, Gold, Black
Model: Red-0704020-028
 Green-0704020-029
 Gold-0704020-030



0202005

Wax loop

Used to remove earwax or dirt in sound outlet

Package: 200pcs/pack Handle Material: ABS

PVC tubing special for IEM



0703046
0703003

Material: PVC
Dimension: 2.5x1.5mm—0703046
 2.0x1.0mm—0703003

Knowles Damper for IEM

100% imported from Knowles, quality assurance IEM.
Several models can be chosen.



Model Variations					
Color	Nom.Acoustic Resistance(Ω)	2.08mm Plug	1.78mm Screen	1.37mm Screen	1.12mm Screen
Grey	330	BF-1999	BF-3034	---	---
White	680	BF-1859	BF-1997	BF-7083	BF-1988
Brown	1000	BF-1860	BF-3035	BF-3292	BF-3320
Green	1500	BF-1861	BF-3038	BF-3039	BF-1991
Red	2200	BF-1921	BF-3037	BF-9778	BF-1995
Orange	3300	BF-1922	BF-3036	---	---
Yellow	4700	BF-1923	---	---	BF-3163



0202011

Brush with wax loop

Used to clean wax or dirty on hearing aids and earmolds

5bunches
Package: 200pcs/pack
Handle Material: ABS
Brush Material: Nylon



0205001

Cleaning cloth

Soft, efficient to wipe earmolds and IEM
Color of cloth will be shipped at random

Dimension: 140x145mm
Material: 80% terylene+20% chinlon



0104030-015

IEM case

Perfect storage for IEM

Material: Aluminum alloy
Color: Black
Dimension: φ80x30mm



0104031

IEM case

Small case for storing IEM, with invisible zipper

Material: EVA
Dimension: φ80x35 mm



0104003-007
0104003-008
0104003-009

IEM case

Perfect storage for IEM

Material: Aluminum alloy
Dimension: φ79x40mm
Color: 0104003-007-Black
0104003-008-Grey
0104003-009-Silver



0104018

IEM case

A small case for storing IEM

Material: Mutispandex
Dimension: 96x70x33 (H)mm



IEM case-S

Used for carrying IEM

Dimension: 90x70x41mm
Color: 0104030-001 Orange
0104030-002 Blue
0104030-003 Transparent
0104030-011 Solid Black



0104022

IEM case

A big case for storing IEM
Two-way zippered closure, easy and secure to use

Material: Oxford
Dimension: 115x85x50 (H)mm



IEM case-M

Dimension: 133x98x40mm
Model: 0104030-004 Transparent Black
0104030-005 Yellow
0104030-006 Blue
0104030-007 Solid Black



0104030-014

IEM case

Perfect storage for IEM

Material: Aluminum alloy
Color: Black
Dimension: φ80x35mm



IEM case-L

Dimension: 140x100x57mm
Model: 0104030-008 Solid Black
0104030-009 Orange
0104030-010 Transparent Black



3D Ear Impression Scanner

Based on structural light technology
Specially developed for customized hearing aid shell& earmold.

- Quickly speed scanning
- High measurement accuracy
- Modular construction
- Efficient and convenient operation

Specifications

Technology	LED raster scanning
Scan time	Around 30s
Number of impression per scanning procedure	1pc
Resolution	Scan data 2 cameras, 130Mega pixels
Accuracy of measuring point	< 12um
Output	STL/ASC/OBJ
Language	Chinese&English
Certificate	CE, FCC
Configuration requirements	Operation system: Windows 7 or Windows 10,64bit CPU: Intel Core i5 or more Graphic card: Independent, NVIDIA GTX640 or GTX750 or higher Memory: 8GB or higher Connect: One USB 3.0/USB 2.0

3D PRINTING SOLUTIONS

Professional and complete solutions including 3D scanner, 3D printer, 3D software, 3D material





DS



MONO



DUO



DUO Scan Touch



Free Stage

3D Impression Scanner

Made In Germany, Fast Scan, High Accuracy

- Without license fees when scanning
- Without any maintenance charge of durable quality
- Removable model holder for comfortable positioning
- Open file formats by STL, ASCII and MSH with any CAD software

Specifications

	DS	DUO	MONO	DUO Scan Touch	Free Stage
Scan time	70-80 seconds	85-95 seconds	90-140 seconds	34-36 seconds	32 seconds
Number of castings per scanning procedure	2	2	1	2	2
Technology	Structured-light scanner with durable high-performance LED	Structured-light scanner with durable high-performance LED	With durable high-performance LED	Structured-light scanner with durable high-performance LED	Structured-light scanner with durable high-performance LED
Weight	10kg	8kg	8kg	8kg	6.5kg
Accuracy of a measuring point	10µm	10µm	10µm	20µm	20µm
Output	STL, ASCII, MSH				
Resolution of scan data	Aadjustable				
Interface	USB				
Power	100-240V AC 50-60Hz, MAX 25W				



3D Printer

Open material system, Fast printing
High-resolution & precision

- High-resolution & precision, layer thickness up to 25-100 µm
- Excellent surface quality
- Faster molding to efficiently complete complex and granular models
- A variety of materials to fit Multi-domain applications
- Compatible with a wide range of brand material supplies

Specifications

Technological Base	DLP
Screen	3.5inch and color IPS touch screen
Light source	385nm LED
Building area	120 * 67.5mm
Output	STL / OBJ / FPP/3MF
Layer thickness	25 ~ 100 µm
Native pixel	+/-50µm
Resolution	1920 * 1080px
Connection	USB.WiFi



3D Printer

Made In Germany
Professional solutions for audio retailers, shops,
labs and producers 30 pieces printed in 35 minutes



- Open material system
- Light source: UV385nm high power LED
- Print thickness: Can adjust to 25-100um
- Output: STL

Specifications

	HA10+	HA20+	HA30+	HA70+
Building area	30 × 60 mm and 60 × 90 mm	133 × 75 mm	151 × 85 mm	232 × 137 mm
Native pixel	+/- 34.5 μm	+/- 34.5 μm	+/- 39 μm	+/- 22.5 μm
Max. part height	90 mm	90 mm	110 mm	150 mm
Resolution	1280 × 720 px	HD 1920 × 1080 px	HD 1920 × 1080 px	4K 3640 × 2160 px XPR
Dimension (W × H × D)	335 × 541 × 349 mm	335 × 541 × 349 mm	480 × 690 × 410 mm	443 × 1593 × 625 mm
Connections	WLAN, TCP/IP, USB	WLAN, TCP/IP, USB	WLAN, TCP/IP, USB	TCP/IP, USB
Control	7" LCD-Display, touch control	7" LCD-Display, touch control	10" touch-screen	10" touch-screen

3D Printer

Model: MAX series , MAX UV
PRO series , PRO 4K



- Open material system
- DLP Technology
- LED light source
- Free and unlimited lifetime technical support
- High-resolution and reliability



Specifications

	MAX UV	PRO 4K 65UV	PRO 4K 80UV
Building area	119 x 67 x 76mm	176.5 x 99 x 200mm	217 x 122 x 200mm
Native pixel	62μm	46μm	6μm
System sizing	260 x 380 x 370mm /16.50Kg	465 x 540 x 1370mm / 140 kg	465 x 540 x 1370mm / 140 kg
File inputs	STL, SLC, STM (Asiga Stomp file format)		
Network Compatibility	Wifi, WirelessDirect, Ethernet		
Software	Asiga Composer software. Lifetime updates included		
Power requirements	100-240VAC, 50/60Hz, 500 Watts (100V - 5Amp Max. 240V - 2.1Amp)		



Secret Ear System

Cyfex is built at Eidgenoessische Technische Hochschule in 2000, in collaboration with leading laboratories. Cyfex precisely develops the features needed so that businesses in dynamic markets of the hearing aid industry are ahead of the game.



Secret Ear Designer®(SED)

Add-on: Milling

The milling tool, available for purchase separately, calculates the optimal position of the earmold in the blank on the basis of the parameters of the milling machine and supports the designers in their work.

The calculated values of the best possible positions and the sound and vent bores are stored in the resulting SPL file that is subsequently transmitted to the CAM software.




Secret Ear Manager®(SEM)

Enables users to effortlessly prepare print jobs.

Among other things orders can be grouped, nested, equipped with support structures and exported as a whole.

Note: A basic version of Secret Ear Manager including the above-mentioned functions is available for free.



Precision is invaluable for laboratories.

The great advantage of generative manufacturing is the high, consistent quality. It is no coincidence that this can be achieved with Cyfex hardware and software.

Many laboratories benefit from the clear concept and streamlined programming of our applications, as well as the optimized interplay of the individual components.



3D Printing Shell material

UV curing premium resin(light spectrum UV 378-388nm) for generative hard ITE shells.

Max construction safety

Suitable for multi-storey print jobs

Accelerated throughput by short light exposure times

Transparent	500g	0404026-001	Transparent	1000g	0404026-015
Beige	500g	0404026-002	Beige	1000g	0404026-016
Intensive red	500g	0404026-003	Intensive red	1000g	0404026-017
Intensive blue	500g	0404026-004	Intensive blue	1000g	0404026-018
Black	500g	0404026-008	Black	1000g	0404026-019
White	500g	0404026-009	White	1000g	0404026-020



3D Printing Hard earmold material

UV curing premium resin(light spectrum UV 378-388nm) for generative manufacturing of hard earmoulds and hearing protection.

Max construction safety

Highest transparency without fillers.

Accelerated throughput by short light exposure times

Transparent	1000g	0404024-010	Blue	1000g	0404024-013
Rose	1000g	0404024-011	Transpantet	1000g	0404027-001
Red	1000g	0404024-012	Transpantet	1000g	0404027-002



3D soft earmold cast material

UV curing resin(light spectrum UV 378-388nm), for generative manufacturing of cast forms, fabrication of silicone earmoulds.

Markedly lower viscosity than conventional material.

Easy to break due to its predefined brittleness, effortless removal of the silicone blank.

Green-transparent	1000g	0404023-003
-------------------	-------	-------------



3D Printing soft earmolds material

UV curing resin (wavelength 385 nm) for generative manufacturing of soft, massive earmoulds

Crystal clear	1000g	0403014-003
---------------	-------	-------------



3D printing cast separator

Separating agent for an effective isolation of cast forms.

Enables an easy removal.

Ready to use dipping solution, transparent, for an efficient use in the lab.

Pack in Aluminum bottle	500ml	0404021-001
-------------------------	-------	-------------

3D Luxaprint shellac

UV curing one component lacquer for coloured sealing of earmoulds made of acrylate.
Colours freely mixable
Easy processing.
Creates a homogeneous scratch-proof, easy to clean surface, reduces cerumen adhesion.



Red	100ml	0404020-001	Black	50ml	0404020-011
Blue	100ml	0404020-002	Red	50ml	0404020-012
Yellow	100ml	0404020-003	Orange	50ml	0404020-013
Orange	100ml	0404020-004	Purple	50ml	0404020-014
Purple	100ml	0404020-005	Green	50ml	0404020-015
Green	100ml	0404020-006	Brown	50ml	0404020-016
Black	100ml	0404020-020	Bue	50ml	0404020-018
Tranparent	300ml	0404020-007	Yellow	50ml	0404020-021

3D Printing hard earmold/shell Material

FotoTec DLP.A is a biocompatible resin with optimal features for the Digital Light Process. Especially for a ~380nm radiation source the material is ideal for manufacturing earmolds and shells.



Optimized initiator system
Low viscosity and minimal adhesive forces
Perfect surface quality
Low viscosity makes cleaning easy

White	500g	0404026-010	Tranparent	1000g	0404024-001
Black	500g	0404026-011	Light orange	1000g	0404024-016
Light blue	1000g	0404026-021	Light purple	1000g	0404024-017
Light red	1000g	0404026-022	Light grey	1000g	0404024-018
Beige	1000g	0404026-005~007	Green	1000g	0404024-019

3D soft earmold cast material

Shortened construction time
Easy cleaning
Easy destroying of cast form after filling and polymerization with material



Light red 1000g 0404023-001

3D printing cast separator

Separating agent for the manufacturing of soft earmolds shell

Valume: 50ml 0404021-002



3D NanoScreen Lacquers

Glazing of all hard earmolds instead of mechanical polishing
Curing is done within 60 seconds.
Brilliant, tissue-friendly surface sealing.

Transparent 250ml 0404020-008
Transparent 20ml 0404020-009



3D FotoScreen Lacquers

The premium lacquer is featured with brilliant surface quality properties on hard earmolds.

Highlights

- ◆ Extremely resistant protective layer for long time
- ◆ Water-clear and extremely glossy
- ◆ Lotus-effect for easy cleaning

Valume: 20ml, 50ml, 250ml
Transparent 0404020-010



0503033-021

3D Dust-free paper

Used to cleaning 3D printed goods.

Size: 21*21cm



0503033-024

Filter nonwoven

Used to filter extra material after printing, specially for illiquid material

Size: 8*15cm



0503033-008

3D Resin case cleaning spray

Used to cleaning 3D printer resin case

Valume:500ml

3D resin material case (Asiga)



0404025-006
0404025-009
0404025-010
0404025-012

Used to storage resin material

Item No. : Used for Pro 1L 0404025-006
10L 0404025-009
5L 0404025-012
Used for Pico 1L 0404025-010

3D resin material case (Rapidshape)



0404025-003
0404025-004
0404025-014

Used to storage resin material

Item No. : 40/60# 0404025-003
30# 0404025-004
20II/30II 0404025-014

Suction machine



0503033-019

Used for extra material resorption

Model: XL01-B
Voltage: AC 220V±22V
Frequency: 50Hz±1Hz
Input: ≤90W
Fuse: RFφ5x20 1.5A/250V
Weight: 4kg
Size: 280 x 190 x 230mm

Material Rotating Device



0503033-025

Special used for color material after storage for long time with material layered.
Easy handle, silent working for earmold lab, hearing aid manufacturer, and 3D lab

Model: XY01-T
Allow weight: 1kg x 2bottle
Speed: 10±3r/min
Size: 230 x 300 x 85mm



0503024-002

3D Connexion

Operate software when 3D printing.
Can realize the design product space selection view in the software,
as well as the far close-range control, local detail viewing



0503033-012

Stainless stell scraper knife

Used to clean 3D printer resin case



0503033-001

3D Manual Centrifuge

Use centrifugal force to remove excess material inside of the printed
earmold or shell, in order to prompt clean and curing.

Material: Plastic
Size: φ230*140mm



0503033-014

3D Material Filter

Can realize the design product space selection view in the software,
as well as the software design product space far close-range control,
local detail viewing

Volume: 1.5L



0503033-013

Reducing valve

Connect nitrogen cylinder and UV curing machine



0503027-009

3D Flash Curing Machine

With two pcs flash lamp for photopolymerization. Special for 3D printing curing for material in 280-580nm.

Brand: NK-Optik
Stroboscopic frequency: 10 flashes/s
Voltage: AC 220V/50HZ
Size: 310 x 310 x 140mm
Capacity: 8pcs ITC or 6pcs ITE
Power: 250W

Model: 3D Flash curing machine 0503027-009
3D Flash curing case 0503029-004 Size: 12*11.9*4.5cm

3D Material & Resin Case Storage

Adjusted temperature, layed display to store earmold material and 3D resin case.

Power: 500W
Temperature: 0-50°C
RH: < 85%
Voltage: 220V±10%
Inner size: 350 x 400 x 500mm
Out size: 540 x 540 x 1040mm

Post-treatment grinder

For wet processing (benchtop machine)

With 1 working container (6 l volume)

- thermally highly cross linking polyurethane working container, inner diameter 230 mm
- concave working container for a high grinding pressure
- centrifugal disc with curved carrier peg for an efficient stream
- precise gap system for thin parts and fine media with a precise distance of 0.1 to 0.3 mm
- exact and robust bearing, low-maintenance
- digital LED-Display for speed control and process time
- adjustment of the process time
- adjustment of the motor speed
- 10 programs for process times
- speed control via frequency inverter
- start/stop for processing flow
- main circuit breaker / emergency stop switch

Core-competitiveness

Soundlink Co., Ltd.

Small MOQ

01

02

Cheap, fast, safe and convenient shipping

Multiple payment terms

03

04

Offer OEM / ODM service

Core values: sincerity seriousness, honesty and thanksgiving

05

06

We regard quality as our life

Manage system is according to ISO 9001:2015

07